

Product datasheet for **SC325687**

ACBD4 (NM_001135707) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ACBD4 (NM_001135707) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACBD4
Synonyms:	HMFT0700
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325687 representing NM_001135707. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCACCGAGAAAGAAAGCCAGAGCCCGACTGCCAGAAACAGTTCCAGGCTGCAGTGAGCGTCATC
CAGAACCTGCCAAGAACGGTCTTACCGCCCTCCTATGAAGAGATGCTGCGATTCTACAGTTACTAC
AAGCAGGCCACCATGGGGCCCTGCCTGGTCCCCGGCCGGGTTCTGGGACCCATTGGACGATATAAG
TGGGACGCCTGGAACAGTCTGGCAAGATGAGCAGGGAGGAGGCCATGTCTGCCTACACTGAAATG
AAACTGGTGGCACAGAAGGTGATCGACACAGTGCCCTGGGTGAGGTGGCAGAGGACATGTTTGGTTAC
TTCGAGCCCTGTACCAGGTGATCCCTGACATGCCGAGGCCCCAGAGACCTTCCGAGAAGGGTCACA
GGTTGAAAGAGCAGGTTGTGAATGGAGATGTTGGGGCTGTTTCAGAGCCTCCCTGCCTCCCAAGGAA
CCGGCACCCCAAGCCAGAGTCCCATTCACCCAGGGACCTGGACTCCGAGGTTTTCTGTGATCCCTG
GAGCAGCTGGAGCCTGAGCTGAGCAGCGGGCAGCATCTGGAGGAAAGCGTGATCCAGGAACAGCCCCG
TGCCCCCACAAGAAAGAGGGGTTGCGGGGCAGCCCGCCGGGGCCCCAGGAGTTGGACGTGTGGCTGC
TGGGGACAGTTTCGAGCACTACAGGAGAGCATGCAGGAGGTGCAGGCGAGGGTGCAGAGCCTGGAGAGCA
TGCCCCGGCCCTGAGCAGTCTCACTCTGTCGGTCAGACTGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001135707

Insert Size: 807 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135707.2](#)

RefSeq Size: 2739 bp

RefSeq ORF: 807 bp

Locus ID: 79777

UniProt ID: [Q8NC06](#)

Cytogenetics: 17q21.31

Protein Families: Transmembrane

MW: 30.3 kDa

Gene Summary: This gene encodes a member of the acyl-coenzyme A binding domain containing protein family. All family members contain the conserved acyl-Coenzyme A binding domain, which binds acyl-CoA thiol esters. They are thought to play roles in acyl-CoA dependent lipid metabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]