

Product datasheet for SC325701

RRNAD1 (NM_001142560) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RRNAD1 (NM_001142560) Human Untagged Clone
Tag:	Tag Free
Symbol:	RRNAD1
Synonyms:	C1orf66; CGI-41; METTL25B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325701 representing NM_001142560. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGGGCATCTCCGCCGAGGCCTCTCTCATGAGGGGAGGAAGCAGCTAGCTGTTAACCTCACCCGT
GTCCTGGCACTCTACCGTTCCATCTTGGATGCCTACATCATCGAATTTTTCACAGACAACCTATGGGAC
ACACTCCCTTGCTCATGGCAGGAAGCATTGGATGGACTGAAACCACCACAGCTGGCCACAATGCTGCTG
GGGATGCCTGGGAAGGGGAGGTCGTACAGGTACAGGTAGTGTGGCCACTCACCTGCTGGCCCTGAAG
TCCACGGCGTGTGCCCTGGCCTTTACCCGGATGCCTGGCTTTAGACCCCTCAGAATTCCTGGAGAAC
CCCAGCCAGAGCTCCCGACTAACAGCTCCATTCCGAAACATGTCAGGCCCAAGAAGCAGCATGAGATC
CGGAGGCTGGGAGAGTTGGTGAAGAAGCTGAGTGATTTACAGGCTGCACCCAGGTTGTAGACGTGGGC
TCAGGCCAGGGCCATCTCTCCCGCTTATGGCTCTTGGCCTGGGTTGATGGTGAAGAGCATCGAAGGG
GATCAGAGACTGGTGGAGAGAGCCAGCGCCTGGACCAGGAGCTTCTGCAGGCTCTGGAGAAAGAGGAG
AAGAGGAACCCGAGATATGTGCAGCGGGGCTACAGCGAGTGGGGCTAGATCCCCAGCTGCCACTGAA
TCTGGCTGCCCTTCAGGCCACGTGGCCAGGAGAACCGTGTGGTGGCCTTCTTCAGCCTGGCTCTACT
GCTTGCCCACTGGTGGAGACGCTTATTCTACTGGACCGCTGCTGTACCTTCAGGAACAGGGTTTCCA
TGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001142560
Insert Size:	834 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142560.1</u>
RefSeq Size:	1798 bp
RefSeq ORF:	834 bp
Locus ID:	51093
UniProt ID:	<u>Q96FB5</u>
Cytogenetics:	1q23.1
MW:	29.7 kDa