

## Product datasheet for SC325752

### MUC15 (NM\_001135092) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	MUC15 (NM_001135092) Human Untagged Clone
Tag:	Tag Free
Symbol:	MUC15
Synonyms:	MUC-15; PAS3; PASIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325752 representing NM_001135092. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCATAATACAATCTATTCTTGCCACATCAAGGGATTGTTATTCCTTTAAAAAAAACCAATACCA
AAGAAGCCTACAATGTTGGCCTTAGCCAAAATTCTGTTGATTTCAACGTTGTTTTATTCATTCTATCG
GGGAGCCATGGAAAAGAAAATCAAGACATAAACACAACACAGAACATTGCAGAAGTTTTTAAAAAATG
GAAAATAAACCTATTTCTTTGGAAAGTGAAGCAAATTAACCTCAGATAAAGAAAATATAACCACCTCA
AATCTCAAGGCGAGTCATTCCCCTCTTTGAATCTACCAACAACAGCCACGGAATAACAGATTTCTCC
AGTAACTCATCAGCAGAGCATTCTTTGGGCAGTCTAAAACCCACATCTACCATTTCCACAAGCCCTCCC
TTGATCCATAGCTTTGTTTCTAAAGTGCCTTGAATGCACCTATAGCAGATGAAGATCTTTTGCCCATC
TCAGCACATCCCAATGCTACACCTGCTCTGTCTTCAGAAAACCTCACTTGGTCTTTGGTCAATGACACC
GTGAAAACCTCTGATAACAGTTCATTACAGTTAGCATCTCTCTTCAGAACCAACTTCTCCATCTGTG
ACCCCTTGATAGTGGAAACCAAGTGGATGGCTTACCACAAACAGTGATAGCTTCACTGGGTTTACCCCT
TATCAAGAAAAACAACCTACAGCCTACCTTAAAATTCACCAATAATTCAAAACCTTTCCAAATACG
TCAGATCCCCAAAAGTTCTGCGATTAGACAATGCACCGAACCTTATGATGTGAGTTTTGGGAATTCT
CCTATGGATGACATACCTCCACTTCGTACTTCTGTATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001135092

**Insert Size:** 936 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001135092.1](#)

**RefSeq Size:** 3259 bp

**RefSeq ORF:** 936 bp

<b>Locus ID:</b>	143662
<b>UniProt ID:</b>	<a href="#">Q8N387</a>
<b>Cytogenetics:</b>	11p14.2
<b>Protein Families:</b>	Secreted Protein, Transmembrane
<b>MW:</b>	33.8 kDa
<b>Gene Summary:</b>	May play a role in the cell adhesion to the extracellular matrix.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. The resulting isoform (c) is shorter than isoform a.