

Product datasheet for **SC325760**

Calumenin (CALU) (NM_001130674) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calumenin (CALU) (NM_001130674) Human Untagged Clone
Tag:	Tag Free
Symbol:	CALU
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325760 representing NM_001130674. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCTGCGACAGTTTCTTATGTGCCTGTCCCTGTGCACAGCCTTTGCCTTGAGCAAACCCACAGAA
AAGAAGGACCGTGTACATCATGAGCCTCAGCTCAGTGACAAGGTTCACAATGATGCTCAGAGTTTGTAT
TATGACCATGATGCCTTCTTGGGTGCTGAAGAAGCAAAGACCTTTGATCAGCTGACACCAGAAGAGAGC
AAGGAAAGGCTTGAATGATTGTAGATAAAATAGACGCGGATAAAGATGGGTTTGTGACGGAGGGGGAG
CTGAAATCCTGGATTAAGCACGCCCAGAAGAAATACATATATGACAATGTTGAAAACCAATGGCAGGAG
TTTGATATGAATCAAGACGGCTTAATCTCCTGGGATGAGTACAGAAACGTGACTTATGGCACTTACCTG
GATGATCCAGATCCTGATGATGGATTTAACTATAAACAGATGATGGTTAGAGATGAGCGGAGGTTTAAA
ATGGCAGACAAGGATGGAGACCTCATTGCCACCAAGGAGGAGTTCACAGCTTTCCTGCACCCTGAGGAG
TATGACTACATGAAAGATATAGTAGTACAGGAAACAATGGAAGATATAGATAAGAATGCTGATGGTTTC
ATTGATCTAGAAGAGTATATTGGTGACATGTACAGCCATGATGGGAATACTGATGAGCCAGAATGGGTA
AAGACAGAGCGAGAGCAGTTTGTGAGTTTCGGGATAAGAACCCTGATGGGAAGATGGACAAGGAAGAG
ACCAAAGACTGGATCCTTCCCTCAGACTATGATCATGCAGAGGCAGAAGCCAGGCACCTGGTCTATGAA
TCAGACAAAACAAGGATGGCAAGCTTACCAAGGAGGAGATCGTTGACAAGTATGACTTATTTGTTGGC
AGCCAGGCCACAGATTTTGGGGAGGCCTTAGTACGGCATGATGAGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001130674

Insert Size: 948 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130674.2](#)

RefSeq Size: 5366 bp

RefSeq ORF: 948 bp

Locus ID: 813

UniProt ID: [O43852](#)

Cytogenetics: 7q32.1

Protein Families: Secreted Protein

MW: 37.1 kDa

Gene Summary: The product of this gene is a calcium-binding protein localized in the endoplasmic reticulum (ER) and it is involved in such ER functions as protein folding and sorting. This protein belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Cab45 and the product of this gene. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

Transcript Variant: This variant (2) has an alternate in-frame exon, compared to variant 1. The resulting isoform (b) is of the same length but has a different internal segment, compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.