

Product datasheet for **SC325776**

BDNF (NM_001143810) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BDNF (NM_001143810) Human Untagged Clone
Tag:	Tag Free
Symbol:	BDNF
Synonyms:	ANON2; BULN2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325776 representing NM_001143810. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGTGGAGCCACCAGTTTTCTCCATGAGTGCACAAGGTTAATCCTTGTTACTACTCAGAATGCTGAG
TTTCTACAGAAAGGGTTGCAGGTCCACACATGTTTTGGCGTCTACCCACACGCTTCTGTATGGCATGAC
TGTGCATCCAGAGAAGGGCTGTGCTGTGTACCTCCAGTTCAGTGAATTTAACAACTGATCCCT
GAAAATGGTTTCATAAAGTTCACCAGGTGAGAAGAGTGATGACCATCCTTTTCTTACTATGGTTATT
TCATACTTTGGTTGCATGAAGGCTGCCCCATGAAAGAAGCAAACATCCGAGGACAAGGTGGCTTGGCC
TACCCAGGTGTGGGACCCATGGGACTCTGGAGAGCGTGAATGGGCCAAGGCAGGTTCAAGAGGCTTG
ACATCATTGGCTGACACTTTTGAACACGTGATAGAAGAGCTGTTGGATGAGGACCAGAAAGTTCGGCCC
AATGAAGAAAACAATAAGGACGCAGACTTGTACACGTCCAGGGTATGCTCAGTAGTCAAGTGCCTTTG
GAGCCTCCTCTTCTTTCTGCTGGAGGAATACAAAAATTACCTAGATGCTGCAAACATGTCATGAGG
GTCCGGCGCCACTCTGACCCTGCCCGCCGAGGGGAGCTGAGCGTGTGTGACAGTATTAGTGAGTGGGTA
ACGGCGGCAGACAAAAGACTGCAGTGGACATGTGGGCGGGACGGTCACAGTCTTGAAAAGGTCCCT
GTATCAAAGGCCAACTGAAGCAATACTTCTACGAGACCAAGTGAATCCCATGGGTTACACAAAAGAA
GGCTGCAGGGGCATAGACAAAAGCATTGGAATCCCAGTGCCGAACACCCAGTCGTACGTGCGGGCC
CTTACCATGGATAGCAAAAAGAGAATTGGCTGGCGATTATAAGGATAGACACTTCTTGTGTATGTACA
TTGACCATTAAGGGGAAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001143810

Insert Size: 990 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001143810.1](#)

RefSeq Size: 4057 bp

RefSeq ORF: 990 bp

Locus ID: 627

UniProt ID:	<u>P23560</u>
Cytogenetics:	11p14.1
Protein Families:	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Transmembrane
Protein Pathways:	Huntington's disease, MAPK signaling pathway, Neurotrophin signaling pathway
MW:	37.1 kDa
Gene Summary:	<p>This gene encodes a member of the nerve growth factor family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature protein. Binding of this protein to its cognate receptor promotes neuronal survival in the adult brain. Expression of this gene is reduced in Alzheimer's, Parkinson's, and Huntington's disease patients. This gene may play a role in the regulation of the stress response and in the biology of mood disorders. [provided by RefSeq, Nov 2015]</p> <p>Transcript Variant: This variant (18), also known as V-VIII, represents use of an alternate promoter and differs in the 5' UTR and 5' coding region, compared to variant 1. The resulting isoform (e) contains a longer N-terminus and lacks a predicted signal peptide compared to isoform a.</p>