

Product datasheet for **SC325785**

BSDC1 (NM_001143890) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: BSDC1 (NM_001143890) Human Untagged Clone
Tag: Tag Free
Symbol: BSDC1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC325785 representing NM_001143890.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAAGGGGAGGACGTGGGATGGTGGCGGAGCTGGCTGCAGCAGAGCTACCAAGCAGTCAAAGAG
AAGGCTCGCCTCTATAGCCTGCAGTCGGACCCAGCAACCTACTGTAATGAACCAGATGGGCCCCGGAA
TTGTTTGACGCCTGGCTTTCCAGTTCTGCTTGGAGGAGAAGAAGGGGGAGATCTCAGAGCTCCTTGTA
GGCAGCCCCTCCATCCGGGCCCTACACCAAGATGGTTCCAGCAGCTGTTTCCATTTCAGAATTCTGG
CATCGGATTTCTATAAAGTCCATCAGTTAGAGCAGGAGCAGGCCCGGAGGGACGCCCTGAAGCAGCGG
GCGGAACAGAGCATCTCTGAAGAGCCCGCTGGGAGGAGGGAAGAGGAGCTCATGGGCATTTACCC
ATATCTCCAAAAGAGGCAAAGGTTCTGTGGCCAAAATTTCTACATTCCCTGAAGGAGAACCTGGCCCC
CAGAGCCCCTGTGAAGAGAATCTGGTGACTTCAGTTGAGCCCCAGCAGAGGTGACTCCATCAGAGAGC
AGTGAGAGCATCTCCCTCGTGACACAGATCGCCAACCCGGCCACTGCACCTGAGGCACGAGTGCTACCC
AAGGACCTGTCCAAAAGCTGCTAGAGGCATCCTTGGAGGAACAGGGCCTGGCTGTGGATGTGGGTGAG
ACTGGACCCTCACCCCTATTCACTCCAAGCCCCTAACGCCTGCTGGCCACACCGCGGCCAGAGCCC
AGGCTCCAGCCAGAGTAGAGACTCTGAGGGAGGAGGCCCCACAGACTTACGGGTGTTGAGCTGAAC
TCGGATAGTGGGAAGTCTACACCCTCCAACAATGGAAAGAAAGGCTCAAGCACGGACATCAGTGAGGAC
TGGGAGAAAGACTTTGACTTGGACATGACTGAAGAGGAGGTGCAGATGGCACTTTCAAAGTGGATGCC
TCCGGGGAGCTGGAAGATGTAGAGTGGGAGGACTGGGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
```

Restriction Sites: Sgfl-MluI
ACCN: NM_001143890
Insert Size: 1008 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001143890.1</u>
RefSeq Size:	3044 bp
RefSeq ORF:	1008 bp
Locus ID:	55108
UniProt ID:	<u>Q9NW68</u>
Cytogenetics:	1p35.1
MW:	37.2 kDa