

Product datasheet for SC325792

PLCXD2 (NM_001134478) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLCXD2 (NM_001134478) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLCXD2
Synonyms:	FLJ31579; phosphatidylinositol-specific phospholipase C, X domain containing 2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325792 representing NM_001134478. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTAGCAGTTAGGAAGGCCAGGAGGAACTCAGGATGGGGACCATCTGCTCCCCAACCCAGCGGG
ACAAAGACATCATCGGAGGTCTGCAATGCCGACTGGATGGCCTCGCTCCCCCTCACCTCCACAACCTC
CCCCTTTCCAATCTGGCAATCCCAGGCTCACATGATTCACTCAGCTACTGGGTGGATGAAAAGTCCCCA
GTGGGGCCTGACCAAAACCAAGCTATCAAACGCCTCGCCAGGATCTCCTTGGTGAAGAAGCTAATGAAG
AAGTGGTCTGTGACTCAGAACCTGACATTCGAGAACAGCTGGAAGCTGGGATCCGCTACTTTGACCTG
CGTGTGTCTTCCAAACCAGGGGATGCCGACCAGGAGATCTACTTCATCCATGGGCTTTTGGCATCAAG
GTCTGGGATGGGCTGATGGAAATTGACTCGTTTCTTACACAGCACCCCCAGGAGATTATCTTCTGGAT
TTCAACCACTTCTATGCCATGGATGAGACCCATCACAATGCCTGGTTCTGCGGATCCAGGAGGCCCTT
GGAAACAAGCTGTGCCAGCCTGCAGTGTGGAAAGTTGACGCTGCGAATCTGTGGGAGAAGAAGTGC
CAGGTTCTTATTTCTACCACTGTCCCTTCTACAAGCAGTACCCCTTCTGTGGCCAGGAAAGAAGATT
CCAGCGCCCTGGGCAACACCACAAGTGTGCGCAAACTAATCCTTCTTGGAGACCACTCTGAGTGAG
CGGGCCTCAGGGGCTCCTTCCATGTCTCCAAGCGATCCTCACCCCAAGAGTGAAGACCATTGCCCGG
GGCTTGGTTGGGGCCTCAAGAACACGCTGGTTCATAGGAATCTTCTGCCATCCTGGACTGGGTGAAA
ACTCAGAAGCCTGGAGCCATGGGTGTCAACATCATCATCTGACTTCGTGGACCTGGTGGACTTTGCT
GCGACTGTATCAAGTTGAATGACCTCCTACAGGAGGACACAGCTCTGGCTAAATGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001134478
Insert Size:	1026 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134478.1</u>
RefSeq Size:	7726 bp
RefSeq ORF:	1026 bp
Locus ID:	257068
Cytogenetics:	3q13.2
MW:	38.6 kDa