

Product datasheet for **SC325799**

MTHFD2 (MTHFD2L) (NM_001144978) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MTHFD2 (MTHFD2L) (NM_001144978) Human Untagged Clone
Tag:	Tag Free
Symbol:	MTHFD2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325799 representing NM_001144978. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACGGTGCCGGTCCGCGGCTTCTCGTGCTCCGCGGCCGCCTTGGCCGAGCGCCGGCGTTGGGCAGA
AGCACAGCACCCCTCCGTAAGGGCACCGGAGAGCCCGGGAGTGCCTCCGGGGCTTTCCGAGCAGCGGT
GTGAGACATGAAGCCATTATTATATCAGGAACCGAAATGGCCAAGCATATCCAGAAAGAAATACAGCGA
GGTGTGGAATCATGGGTTCCCTTGGAAACAGAAGACCTCACCTCAGTATAATTTTAGTGGGAGATAAC
CCAGCAAGCCATACATATGTCAGGAATAAGATAAGAGCTGCCTCTGCTGTAGGTATTTGTAGTGAGCTC
ATTCTAAAACCTAAGGATGTTTCTCAGGAAGAACTTTGGACGTAACGATCAATTGAATATGGACCCA
AGAGTCAGCGGTATATTAGTTCAGTTACCACTACCAGACCAGTTGATGAGCGAACAATATGCAATGGA
ATTGCCCCAGAAAAAGATGTAGATGGATTTTCAATTATCAATATTGGAAGATTGTCCTTATCAGCAT
TCTCTCATACCTGCCACTGCCAGTGCTGTTTGGGAAATAATAAAAAGAACAGGAATTCAAACATTTGGA
AAAAATGTGGTTGTGGCTGGAAGATCCAAGAACGTAGGGATGCCTATTGCCATGCTTTTACACACTGAT
GGAGAGCATGAACGCCAGGAGGTGATGCAACTGTGACAATAGCTCACAGATACACCCCAAGAGCAA
CTGAAGATTCATACGCAGCTGGCAGATATTATCATAGTTGCTGCAGGTATCCAAAGTTGATTACGCTC
GATATGGTTAAAGAAGGTGCTGCTGTAATTGATGTGGGTATCAACTATGTCCACGATCCAGTGACAGGA
AAGACAAAATTAGTTGGAGATGTGGACTTCGAAGCTGTTAAAAAGAAAGCTGGCTTTATCACTCCAGTT
CCAGGAGGTGTGGACCCATGACAGTGGCAATGCTTCTGAAGAACACCCCTTCTGGCAGCTAAAAAATC
ATTTACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001144978

Insert Size: 1044 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001144978.1](#)

RefSeq Size: 2372 bp

RefSeq ORF: 1044 bp

Locus ID: 441024
UniProt ID: [Q9H903](#)
Cytogenetics: 4q13.3
MW: 37.3 kDa