

Product datasheet for **SC325809**

TCP10L2 (NM_001145121) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TCP10L2 (NM_001145121) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCP10L2
Synonyms:	bA517H2.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325809 representing NM_001145121. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCTGCTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGCTGGAGGTCAGCTCGAGGCCAGGGAGCCCAAGGAGGGCACCCACCCAGAGGACCCGTGCCGGGA
 GCTGGGGTGGCATGGAGAAGACACCTGCAGCAGCCGAGGTCCCCAGGGAGGACAGCAATGCCGGGGAG
 ATGCCGTCATTACAGCAGCAGATCACCAGCCTCCACCAAGAGCTCGGGAGACAGCAGTCGCTGTGGGCT
 GATATTCACAGAAAACCTCAGAGTCATATGGATGCCTTGAGGAAGCAGAACCGGGAGCTCCGAGAAGAG
 CTGAGAGGCCTGCAGCGCAGCAGTGGGAAGCCGGGAAGAAACCCGAGCGTCCCACACGCGGGGCGA
 GAATCACACACTCTGGCATTGGAACCTGCTTTTGAAAAATTTACCTCTGTGAGCTGATGAAGAGACA
 ACGCCCAAATACGCTGGCCGCAAGAGTCAGAGTGCGACTCTCCTGGGACAAAGATGGTCATCTAACCAC
 TTAGCTCTCTCAAAGCCAATGAGTTTAAAGACAGAAAGAATTAACCTCGGGGAAAAACACCACACAGGAA
 GATAGAGAGAAAAAGTCTCCCGGGAGACGTCAAGACAGAAAGTCCAGCACCCACTGGAAGGCCGACTCCC
 GGTGCAGAAAGACGGGAGGTGTCTGAAGACGGAAGATTATGCATCCATCTTCCCGAAGTCTGCAAAAC
 TCCGGTGGCAGAAATCACCAGTGCAGGCTTCCAGGCCGCCACGCTGCAGGAGCAGACGGCAGCAGCC
 GGAGTTGCTGACAGAAGCTCATCAGTCTTAGGGAGTTCTGAAGGCGGATTCTCAGCCGCGTTCAAGCT
 GAAGAGTTCGCCAGTTCTTCCCCAGACAGTGCAGAGCGGCAGAAATCTTCTGTAAACCCACCATCTTCC
 CTGGAAATAGCACAGGCCATGGACACTAAGATGAAAAAGGAGGAAGTACAGGAAGAAAAACGACATCCA
 AATGGCAAGGCAGATGACTGCCGGCGGTCCGGTTTCCCGAGTGAATTTCCAGGAGCTCTGCATGCCGCT
 CCCTCCAGACAGGACATGGGGCCCTGA
 ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001145121


[View online »](#)

Insert Size:	1062 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145121.1</u>
RefSeq Size:	2185 bp
RefSeq ORF:	1062 bp
Locus ID:	401285
Cytogenetics:	6q27
MW:	38.4 kDa