

Product datasheet for **SC325847**

SCRN2 (NM_001145023) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	SCRN2
Synonyms:	Ses2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325847 representing NM_001145023. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGTCGTCGAGCCCTGACTCCCATGTTCTGCGACTGCTTTGTCTCCGTGCCCGGCCCTCAGCC
ATCCCGGCTGTGATCTTTGCCAAGAACTCGGACCGACCCCGGGACGAGGTGCAGGAGGTGGTGTGTTGTC
CCCGCAGGCACTCACACTCCTGGGAGCCGGCTCCAGTGCACCTACATTGAAGTGAACAGGTGTCTGAAG
ACGCACGCTGTGATTCTGAGCCGCTCTTCTGGCTATGGGGGGCTGAGATGGGCGCCAACGAGCATGGT
GTCTGCATTGGCAACGAGGCTGTGTGGACGAAGGAGCCAGTTGGGAGGGGGAAGCCCTGCTGGGCATG
GACCTACTCAGGCTGGCTTTGGAACGGAGCAGCTCTGCCAGGAGGCCTTGATGTGATCACAGGGTTA
CTGGAGCACTATGGGAGGGGGGCAACTGCCTGGAGGATGCTGCGCCATTCTCTACCATAGCACCTTC
CTGCTGGCTGACCGCACTGAGGCGTGGGTGCTGGAGACAGCTGGGAGGCTCTGGGCTGCACAGAGGATC
CAGGAGGGGGCCCAACATCTCCAACAGCTGAGCATTGGCACGGACATCTCGGCCCAACACCCGGAG
CTGCGGACTCATGCCAGGCCAAGGGCTGGTGGGATGGGAGGGTGCCTTTGACTTTGCTCAGATCTTC
TCCCTGACCCAGCAGCCTGTGCGCATGGAGGCTGCCAAGGCCCGCTTCCAGGCAGGGCGGGAGCTGCTG
CGGCAACGGCAAGGGGGCATCACGGCAGAGGTGATGATGGGCATCCTCAGAGACAAGGAGAGTGGTATC
TGATGGACTCGGAGGCTTTTCGACCACGGCCAGCATGGTGTCTGTCTGCCCCAGGATCCACGCAG
CCCTGCGTGCACTTTCTTACCGCCACGCCAGACCATCCAGGTCTGTGTTCAAACCTTTTCATCTTCGGG
ATGGGGGTGGCCAGGCCCGCCAGGTGCTGTCCCCCACTTTTGGAGCACAAGACCTGTTTCGACCCCTG
CCCCGATTCCAGACTCAGGTAGATCGTCGGCATACCTCTACCGTGGACACCAGGCAGCCCTGGGGCTG
ATGGAGAGAGATCAGGTATCCCCAGGGAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



ACCN:	NM_001145023
Insert Size:	1137 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001145023.1</u>
RefSeq Size:	1849 bp
RefSeq ORF:	1137 bp
Locus ID:	90507
UniProt ID:	<u>Q96FV2</u>
Cytogenetics:	17q21.32
MW:	41.3 kDa