

Product datasheet for **SC325849**

ZNF669 (NM_001142572) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF669 (NM_001142572) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF669
Vector: pCMV6 series
Fully Sequenced ORF: >NCBI ORF sequence for NM_001142572, the custom clone sequence may differ by one or more nucleotides

```

ATGGACTCGGTGGCTTTTGAGGATGTGGCTGTGAACTTTACCCAGGAGGAATGGGCTTTG
CTAGATTCTTCTCAGAAGAATCTCTACAGAGAAGTGATGCAGGAAACCTGCAGGAACCTG
GCTTCTGTAGGAAGCCAATGGAAAGACCAGAATATTGAAGATCACTTCGAAAAACCTGGG
AAAGATATAAGAAATCATATCGTACAGAGACTGTGTGAAAGTAAAGAAGATGGTCAGTAT
GGAGAAGTTGTCAGCCAAATCCAAATCTTGATCTGAACGAGAACATTTCTACTGGATTA
AAACCATGTGAATGCAGTATTTGTGGAAAAGTCTTTGTACGTCATTCCCTCCTTAATAGG
CATATCCTAGCTCACTCAGGATACAAACCATATGGAGAGAAGCAATATAAATGTGAACAG
TGTGGGAAATTCTTCGTTTCTGTTCCAGGTGTTAGAAGACACATGATAATGCACAGTGGA
AATCCAGCTTATAAATGTACGATATGTGGGAAAGCTTTTTATTTTCTCAATTCAGTTGAA
AGACATCAGAGAAGTACACAGGAGAGAAAAACCTATAAATGTAAACAATGTGGTAAAGCA
TTCAGTGTTCGGTCTTGTCTAATACATGAACGAACTCACACTGGAGAGAAACCTAC
GAATGTAAGGAATGTGGGAAAACATTCAGATTTTCTTGTCTTTTAAGACGCATGAAAGG
ACTCACACTGGAGAAAGACCTATAAATGTACCAAATGTGATAAAGCCTTCAGCTGTTCC
ACTTCCCTTCGTTACCATGGAAGCATTCTACTGGAGAGAGACCCTATGAGTGTAACAA
TGTGGCAAAGCCTTTAGTCGTTTGAGTTCCTTTGTAACCATAGAAGTACTCATACCGGA
GAGAAACCTATGAATGTAAACAATGTGATCAAGCCTTCAGTCGCCTCAGTTCCTTCAC
CTCCACGAAAGAAATTCATACTGGAGAAAAACCTATGAATGTAAGAAATGCGGTAAAGCC
TACACTCGTTCAGTCACCTTACTCGCCATGAAAGAAGTCATGATATAGAGGCTGGGTGT
AGTGACTCAGCTATAATCCCAGCACTTTGGGAGGCCAAGGCGTGTGGATTGCT
  
```

Restriction Sites: Please inquire
ACCN: NM_001142572
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142572.1</u> , <u>NP_001136044.1</u>
RefSeq Size:	1740 bp
RefSeq ORF:	1137 bp
Locus ID:	79862
UniProt ID:	<u>Q96BR6</u>
Cytogenetics:	1q44
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 5' coding region, compared to variant 1, which results in a protein (isoform 2) that lacks a segment of the N-terminus, compared to isoform 1.