

Product datasheet for **SC325871**

BIGM103 (SLC39A8) (NM_001135148) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: BIGM103 (SLC39A8) (NM_001135148) Human Untagged Clone
Tag: Tag Free
Symbol: SLC39A8
Synonyms: BIGM103; CDG2N; LZT-Hs6; PP3105; ZIP8
Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_001135148, the custom clone sequence may differ by one or more nucleotides

```
ATGCACCAACACGCTCAGTGTAACTGCTGAAGAGATCTTTCCCTTCATGGCTTTTCA
AATGCTACCCAAATAACCGCTCCAAATCTCTGTCTGTCCAGCAGTCTTACAGCAA
TTGAACTTTCACCCATGTGAGGATCGGCCCAAGCACAAAACAAGACCAAGTCATTAGAA
GTTTGGGGATATGGATTCTGTGAGTACGATTATTAATCTGGCATCTCTCCTCGGATTG
ATTTTGACTCCACTGATAAAGAAATCTTATTTCCCAAAGATTTTGACCTTTTTTGTGGG
CTGGCTATTGGGACTCTTTTTTCAAATGCAATTTTCCAATTATCCAGAGGCATTGGA
TTTGATCCCAAAGTCGACAGTTATGTTGAGAAGGCAGTTGCTGTGTTGGTGGATTTAC
CTACTTTTCTTTTTGAAAGAATGCTAAAGATGTTATTAAGACATATGGTCAGAATGGT
CATACCCACTTTGGAAATGATACTTTGGTCTCAAGAAAAAAGTCAACCTAAAGCA
TTACCTGCCATCAATGGTGTGACATGCTATGCAAATCCTGCTGTACAGAAAGCTAATGGA
CATATCCATTTTGATAATGTGAGTGTGGTATCTCTACAGGATGGAAAAAAGAGCCAAAGT
TCATGTACCTGTTTGAAGGGGCCAAACTGTGAGAAATAGGGACGATTGCCTGGATGATA
ACGCTCTGCGATGCCCTCCACAATTTTCATCGATGGCCTGGCGATTGGGGCTTCTGCACC
TTGTCTCTCCTTCAGGGACTCAGTACTTCCATAGCAATCCTATGTGAGGAGTTTCCCAC
GAGTTAGGAGACTTTGTGATCTACTCAATGCAGGGATGAGCACTCGACAAGCCTTGCTA
TTCAACTTCCTTTCTGCATGTTCCCTGCTATGTTGGGCTAGCTTTTGGCATTGTTGGTGGC
AACAAATTCGCTCCAAATATTATATTTGCACTTGTGGAGGCATGTTCTCTATATTTCT
CTGGCAGATATGTTCCAGAGATGAATGATATGCTGAGAGAAAAGGTAAGTGAAGAAAA
ACCGATTTACCTTCTTCATGATTCAGAATGCTGGAATGTTAACTGGATTCACAGCCATT
CTACTATTACCTTGTATGCAGGAGAAATCGAATTGGAG
```

Restriction Sites: Please inquire

ACCN: NM_001135148

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001135148.1</u> , <u>NP_001128620.1</u>
RefSeq Size:	2624 bp
RefSeq ORF:	1182 bp
Locus ID:	64116
UniProt ID:	<u>Q9C0K1</u>
Cytogenetics:	4q24
Protein Families:	Transmembrane
Gene Summary:	<p>This gene encodes a member of the SLC39 family of solute-carrier genes, which show structural characteristics of zinc transporters. The encoded protein is glycosylated and found in the plasma membrane and mitochondria, and functions in the cellular import of zinc at the onset of inflammation. It is also thought to be the primary transporter of the toxic cation cadmium, which is found in cigarette smoke. Multiple transcript variants encoding different isoforms have been found for this gene. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (4) has a distinct 5' exon, compared to variant 1. The resulting protein (isoform c) has a shorter and distinct N-terminus, compared to isoform a.</p>