

## Product datasheet for **SC325898**

### ZFP92 (NM\_001136273) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZFP92 (NM\_001136273) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZFP92  
**Synonyms:** ZNF897  
**Vector:** pCMV6 series

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001136273, the custom clone sequence may differ by one or more nucleotides

```
ATGGCAGCCATTCTCCTGACCACGAGACCCAAGGTGCCAGTATCTTTTGAGGATGTGTCC
GTGTACTTCACAAAGACAGAATGGAAGCTTCTGGACCTCAGACAAAAGTCCTCTACAAG
CGGGTGATGCTGGAGAATATAGCCATTTGGTGTCACTGGGATTTTCTTCTCCAAGCCT
CACCTGATCTCCAATTGGAACGAGGGGAAGGACCCTGGGTAGCAGACATCCCCAGAACC
TGGGCCACCGCAGGATTGCACATAGGTGACAGAACACAGAGCAAGACGTCGACTTCAACG
CAGAAGCATTCTGGACGACAACCTCCCCGGGCGCATCCACAAGGTGGCAAGGAGGGGCGAG
GCGGGCAGGTCGTCTGTGCTCCAGAGAGGTGCCAGGGCTTGGGGCAGAGTTCGGCTGCG
GGGCCGCAGGGCCCCAAAGGCGCGGAGAAGCGGTACCTGTGCCAGCAGTGTGGGAAGGCC
TTCAGCCGCGAGCTCCAACCTCATCAAGCACCGCATCATCCACAGTGGCGAGAAGCCTTAC
GCGTGCCCCGAGTGCGGCAAGCTGTTTCGCCGAGCTTCGCGCTCCTGGAGCACCAGCGC
ATCCACAGCGGCGAGAAGCCCTACGCTGCCCGAGTGCAGCAAGACCTTACGCGCAGC
TCCAACCTCATCAAGCACAGGTCATCCACAGCGGCGAGCGGCCCTTCGCTGCGGCGAC
TGCGGCAAAGTGTCCGCGCAGCTTCGCGCTCCTGGAGCACGCGCGCTGCACAGCGGC
GAGCGGCCCTACGCGTGCCAGAGTGCAGCAAGGCCTTACGCCGAGCTCCAACCTCATC
GAGCACCAGCGCACGACCGCGGCGAGAAGCCCTACGCTGCGGCCAGTGCGCCAAGGCC
TTCAAGGGCGTCTCGCAGCTCATCCACCACAGCGCAGCCACAGCGCGAGCGGCCCTTC
GCGTGCCGCGAGTGCGGCAAGGCCTTCCGTGGCCGTTCCGGCCTCAGCCAGCACCAGCGC
GTGCACAGCGGTGAGAAGCCCTACGAGTGCAGCGACTGCGGCAAGGCCTTCGGCCGGCGC
GCCAACCTATTCAAGCACAGGCGAGTGCAGCGGCCAGGCGCCCTGCGAAGGCGGAGACG
GCGCGGAGGCTAGCGGGCCCTGGGAGCACCGGCCCTGGGAGCGCGGTGGCGGCCACCAGC
CCCCGCGGCCGAGCACAGCCGCGAGGCCTTCCAGGCCAGCCGCGCGC
```

**Restriction Sites:** Please inquire

**ACCN:** NM\_001136273

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001136273.1</a> , <a href="#">NP_001129745.1</a>
<b>RefSeq Size:</b>	1251 bp
<b>RefSeq ORF:</b>	1251 bp
<b>Locus ID:</b>	139735
<b>UniProt ID:</b>	<a href="#">A6NM28</a>
<b>Cytogenetics:</b>	Xq28
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]