

Product datasheet for **SC325921**

TMEM130 (NM_001134450) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM130 (NM_001134450) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM130
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001134450, the custom clone sequence may differ by one or more nucleotides

```

ATGGCCCAGGCAGTGTGGTCGCGCCTCGGCCGATCCTCTGGCTTGCCTGCCTCCTGCC
TGGGCCCCCGGCAGGGGTGGCCGAGGCCTGTATGAACTCAATCTCACCACCGATAGCCCT
GCCACCACGGGAGCGGTGGTGACCATCTCGGCCAGCCTGGTGGCAAGGACAACGGCAGC
CTGGCCCTGCCGCTGACGCCACCTCTACCGCTTCCACTGGATCCACACCCCGTGGTG
CTTACTGGCAAGATGGAGAAGGGTCTCAGCTCCACCATCCGTGTGGTGGCCACGTGCC
GGGGAATTCCCGTCTCTGTCTGGTCACTGCCGCTGACTGCTGGATGTGCCAGCCTGTG
GCCAGGGCTTTGTGGTCTCCCATCACAGAGTTCTCGTGGGGACCTTGTGTGCC
CAGAACACTTCCCTACCCTGGCCAGCTCCTATCTCACTAAGACCGTCCGAAAGTCTCC
TTCTCTCCACGACCCGAGCAACTTCTCAAGACCGCCTTGTCTCTACAGCTGGGAC
TTCGGGGACGGGACCCAGATGGTGAAGACTCCGTGGTCTATTATACTATTCCATC
ATCGGGACCTTACCCTGAAGCTCAAAGTGGTGGCGGAGTGGGAAGAGGTGGAGCCGAT
GCCACGAGGGCTGTGAAGCAGAAGACCGGGGACTTCTCCGCTCGCTGAAGCTGCAGGAA
ACCTTCGAGGCATCAAGTGTGGGGCCACCCTAATTCAGACCTTCCAAAAGATGACC
GTGACCTTGAAGTCTCTGGGAGCCCTCCTCTGACTGTGTGCTGGCGTCTCAAGCCTGAG
TGCTCCCGCTGGAGGAAGGGGAGTGCCACCCTGTGTCCGTGGCCAGCACAGCGTACAAC
CTGACCCACACCTTCAGGACCCCTGGGGACTACTGCTTCAGCATCCGGGCCGAGAATATC
ATCAGCAAGACACATCAGTACCACAAGATCCAGGTGTGGCCCTCCAGAATCCAGCCGGCT
GTCTTTGCTTTCCCATGTGCTACACTTATCACTGTGATGTTGGCCTTCATCATGTACATG
ACCCTGCGGAATGCCACTCAGCAAAAGGACATGGTGGAGGTGGCTGATTTTGACTTTTCC
CCCATGTCTGACAAGAACCCGGAGCCACCCTCTGGGGTCAGGTGCTGCTGCCAGATGTGC
TGTGGGCTTTCTTGTGGAGACTCCATCTGAGTACCTGGAAATTGTTCTGTGAGAACCAC
GGGCTGCTCCCGCCCTCTATAAGTCTGTCAAACTTACACCGTG

```

Restriction Sites:	Please inquire
ACCN:	NM_001134450
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001134450.1 , NP_001127922.1
RefSeq Size:	3065 bp
RefSeq ORF:	1308 bp
Locus ID:	222865
UniProt ID:	Q8N3G9
Cytogenetics:	7q22.1
Protein Families:	Transmembrane