

Product datasheet for SC325922

MEF2A (NM_001130927) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MEF2A (NM_001130927) Human Untagged Clone
Tag:	Tag Free
Symbol:	MEF2A
Synonyms:	ADCAD1; mef2; RSRFC4; RSRFC9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325922 representing NM_001130927. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGCGGAAGAAAATACAAATCACACGCATAATGGATGAAAGGAACCGACAGACTTTAAGAAAGAAA
GGCCTTAATGGTTGTGAGAGCCCTGATGCTGACGATTCTTTGAGCACAGTCCACTCTCGGAGGACAGA
TTCAGCAAATAAATGAAGATAGTATTTTATTTTCAAACGAGGCCCTCTGGTCTGCCACCTCAGAAC
TTTTCAATGTCTGTCACAGTTCAGTGACCAGCCCAATGCTTTGCTCTACACTAACCCAGGGAGTTCA
CTGGTGTCCCATCTTTGGCAGCCAGCTCAACGTTAACAGATTCAAGCATGCTCTCTCCACCTCAAACC
ACATTACATAGAAATGTGTCTCCTGGAGCTCCTCAGAGACCACCAAGTACTGGCAATGCAGTGGGATG
TTGAGCACTACAGACCTCACAGTGCCAAATGGAGCTGGAAGCAGTCCAGTGGGGAATGGATTTGTAAC
TCAAGAGCTTCTCCAAATTTGATTGGAGCTACTGGTGAAATAGCTTAGGCAAATCATGCCCTACAAG
TCTCCCCCTCCACCAGGTGGTGGTAATCTTGAAATGAACAGTAGGAAACCAGATCTTCGAGTTGTCATC
CCCCCTCAAGCAAGGGCATGATGCCTCCACTATCGGAGGAAGAGGAATTGGAGTTGAATACCCAGAGG
ATCAGTAGTTCTCAAGCCACTCAACCTCTTGCTACCCAGTCTGTGTGACAACCCCAAGCTTGCCCT
CCGCAAGGACTTGTGACTCAGCAATGCCGACTGCCTACAACACTGATTATCACTGACCAGCGCTGAC
CTGTGAGCCCTTCAAGGCTTCAACTCGCCAGGAATGCTGTGCTGGGACAGGTGTCGGCTGGCAGCAG
CACCACCTAGGACAAGCAGCCCTCAGCTCTTTGTTGCTGGAGGGCAGTTATCTCAGGGTTCCAATTTA
TCCATTAATACCAACCAAAACATCAGCATCAAGTCCGAACCGATTTCACCTCCTCGGGATCGTATGACC
CCATCGGGCTTCCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAG
CCACAACCCCGCAGCCCGACAGGAAATGGGGCGCTCCCCTGTGGACAGTCTGAGCAGCTCT
AGTAGCTCCTATGATGGCAGTGATCGGGAGGATCCACGGGGCGACTTCCATTCTCCAATTGTGCTTGGC
CGACCCCAAACTGAGGACAGAGAAAGCCCTTCTGTAAGCGAATGAGGATGGACGCGTGGGTGACC
TAA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTAA
```



[View online >](#)

Restriction Sites: Sgfl-RsrII

Plasmid Map:



ACCN: NM_001130927

Insert Size: 1314 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130927.2](#)

RefSeq Size: 5368 bp

RefSeq ORF: 1314 bp

Locus ID: 4205

UniProt ID: [Q02078](#)

Cytogenetics: 15q26.3

Protein Families: Transcription Factors

MW: 46.5 kDa

Gene Summary: The protein encoded by this gene is a DNA-binding transcription factor that activates many muscle-specific, growth factor-induced, and stress-induced genes. The encoded protein can act as a homodimer or as a heterodimer and is involved in several cellular processes, including muscle development, neuronal differentiation, cell growth control, and apoptosis. Defects in this gene could be a cause of autosomal dominant coronary artery disease 1 with myocardial infarction (ADCAD1). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2010]

Transcript Variant: This variant (3) lacks an in-frame 5' coding exon and a 5' non-coding exon, compared to transcript variant 6. These differences result in a shorter isoform (3), compared to isoform 5. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.