

Product datasheet for **SC325946**

TRIP15 (COPS2) (NM_001143887) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIP15 (COPS2) (NM_001143887) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIP15
Synonyms:	ALIEN; CSN2; SGN2; TRIP15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325946 representing NM_001143887. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTGACATGGAGGATGATTTTCATGTGCGATGATGAGGAGGACTACGACCTGGAATACTCTGAAGAT
AGTAACTCCGAGCCAAATGTGGATTTGGAAAATCAGTACTATAATTCCAAAGCATTAAAAGAAGATGAC
CCAAAAGCGGCATTAAGCAGTTTCCAAAAGTTTTGAACTTGAAGGTGAAAAGGAGAATGGGGATTT
AAAGCACTGAAACAAATGATTAAGATTAACCTCAAGTTGACAACTTTCCAGAAATGATGAATAGATAT
AAGCAGCTATTGACCTATATTCGGAGTGCAGTCACAAGAAATTTCTGAAAAATCCATTAATTCTATT
CTTGATTATATCTACTTCTAAACAGAATTCTGATTTTTTATGTCAGATGGATTTACTGCAGGAATTC
TATGAAACAACACTGGAAGCTTTGAAAGATGCTAAGAATGATAGACTGTGGTTTAAAGACAAACAAAG
CTTGAAAATTATATTTAGAACGAGAGGAATATGGAAAAGCTTCAAAAAATTTTACGCCAGTTACATCAG
TCGTGCCAGACTGATGATGGAGAAGATGATCTGAAAAAGGTACACAGTTATTAGAAATATATGCTTTG
GAAATTCAAATGTACACAGCACAGAAAAATAACAAAAAAGCTTAAAGCACTCTATGAACAGTCACTTCAC
ATCAAGTCTGCCATCCCTCATCCACTGATTATGGGAGTTATCAGAGAATGTGGTGGTAAAATGCACCTG
AGGGAAGGTGAATTTGAAAAGGCACACACTGATTTTTTTGAAAGCCTTCAAGAATTATGATGAATCTGGA
AGTCCAAGACGAACCACTTGCTTAAAATATTTGGTCTTAGCAAAATGCTTATGAAATCGGGAATAAAT
CCATTTGACTCACAGGAGGCCAAGCCGTACAAAAATGATCCAGAAATTTAGCAATGACGAATTTAGTA
AGTGCCTATCAGAATAATGACATCACTGAATTTGAAAAGATTCTAAAAACAATCAGCAACATCATG
GATGATCCTTTTATAAGAGAACACATTGAAGAGCTTTTGCGAAACATCAGAACACAAAGTCTTATAAAA
TTAATTAAGCCTTACACAAGAATACATATTCCTTTTATTTCTAAGGAGTTAAACATAGATGTAGCTGAT
GTGGAGAGCTTGCTGGTGCAGTGCATATTGGATAACACTATTCATGGCCGAATTGATCAAGTCAACCAA
CTCCTTGAAGTGGATCATCAGAAGAGGGTGGTGCACGATATACTGCACTAGATAAATGGACCAACCAA
CTAAATCTCTCAACCAGGCTGTAGTCAGTAAACTGGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



View online »

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001143887
Insert Size:	1353 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001143887.1</u>
RefSeq Size:	4124 bp
RefSeq ORF:	1353 bp
Locus ID:	9318
UniProt ID:	<u>P61201</u>
Cytogenetics:	15q21.1
Protein Families:	Druggable Genome, Transcription Factors
MW:	52.4 kDa
Gene Summary:	<p>Essential component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. The complex is also involved in phosphorylation of p53/TP53, c-jun/JUN, IκappaBα/NFKBIA, ITPK1 and IRF8/ICSBP, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively. Involved in early stage of neuronal differentiation via its interaction with NIF3L1.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) represents the longer transcript and encodes the longer isoform (2).</p>