

## Product datasheet for **SC325956**

### SGK1 (NM\_001143677) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SGK1 (NM_001143677) Human Untagged Clone
Tag:	Tag Free
Symbol:	SGK1
Synonyms:	SGK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC325956 representing NM\_001143677.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCATCTCAGAGTCCAGCTTATCAGAGGCATGTAGCAGGGAGGCTTATTCCAGCCATAACTGGGCT
CTACCTCCAGCCTCCAGAAGTAATCCCCAACCTGCATATCCTTGGCAACCCGAAGAATGAAAGAAGAA
GCTATAAAACCCCTTTGAAAGCTTTCATGAAGCAGAGGAGGATGGGTCTGAACGACTTTATTGAGAAG
ATTGCCAATAACTCCTATGCATGCAAACACCCTGAAGTTCAGTCCATCTTGAAGATCTCCCAACCTCAG
GAGCCTGAGCTTATGAATGCCAACCTTCTCCTCCACCAAGTCTTCTCAGCAAATCAACCTTGGCCCG
TCGTCCAATCCTCATGCTAAACCATCTGACTTTCACTTCTTGAAGTGATCGAAAGGGCAGTTTTGGA
AAGTTCTTCTAGCAAGACACAAGGCAGAAGAAGTGTCTATGCAGTCAAAGTTTACAGAAGAAAGCA
ATCCTGAAAAAGAAAGAGGAGAAGCATATTATGTCGGAGCGGAATGTTCTGTTGAAGAATGTGAAGCAC
CCTTCTGTTGGCCTTCACTTCTTTCCAGACTGCTGACAAATTGTACTTTGCTAGACTACATT
AATGGTGGAGAGTTGTTCTACCATCTCCAGAGGAACGCTGCTTCTGGAACCACGGGCTCGTTTCTAT
GCTGCTGAAATAGCCAGTGCCTTGGGCTACCTGCATTCACTGAACATCGTTTATAGAGACTTAAACCA
GAGAATATTTTGCTAGATTACACAGGGACACATTGTCCTTACTGACTTCGGACTCTGCAAGGAGAACATT
GAACACAACAGCACAACATCCACCTTCTGTGGCAGCCGGAGTATCTCGCACCTGAGGTGCTTCATAAG
CAGCCTTATGACAGGACTGTGGACTGGTGGTGCCTGGGAGCTGTCTTGTATGAGATGCTGTATGGCCTG
CCGCCTTTTATAGCCGAAACACAGCTGAAATGTACGACAACATTCTGAACAAGCCTCTCCAGCTGAAA
CCAAATATTCAAATCCGCAAGACACCTCCTGGAGGGCTCCTGCAGAAGGACAGGACAAAGCGGCTC
GGGCAAGGATGACTTCATGGAGATTAAGAGTCATGTCTTCTCCTTAATTAAGTGGGATGATCTC
ATTAATAAGAAGATTACTCCCCCTTTAAACCAATGTGAGTGGGCCCAACGACCTACGGCACTTTGAC
CCCGAGTTTACCGAAGAGCCTGTCCCAACTCCATTGGCAAGTCCCTGACAGCGTCCTCGTCACAGCC
AGCGTCAAGGAAGCTGCCGAGGCTTTCCTAGGCTTTCCTATGCGCTCCACGGACTCTTTCCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



**ACCN:** NM\_001143677

**Insert Size:** 1380 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001143677.1</a></u>
<b>RefSeq Size:</b>	2459 bp
<b>RefSeq ORF:</b>	1380 bp
<b>Locus ID:</b>	6446
<b>UniProt ID:</b>	<u><a href="#">O00141</a></u>
<b>Cytogenetics:</b>	6q23.2
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>MW:</b>	52.1 kDa
<b>Gene Summary:</b>	<p>This gene encodes a serine/threonine protein kinase that plays an important role in cellular stress response. This kinase activates certain potassium, sodium, and chloride channels, suggesting an involvement in the regulation of processes such as cell survival, neuronal excitability, and renal sodium excretion. High levels of expression of this gene may contribute to conditions such as hypertension and diabetic nephropathy. Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jan 2009]</p> <p>Transcript Variant: This variant (3) contains an alternative in-frame, 5' terminal exon compared to transcript variant 1, resulting in an isoform (3) with a longer and an unique N-terminus compared to isoform 1.</p>