

Product datasheet for **SC325971**

ARFGAP3 (NM_001142293) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARFGAP3 (NM_001142293) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARFGAP3
Synonyms:	ARFGAP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC325971 representing NM_001142293.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGGACCCAGCAAGCAGGACATCTTGACCATCTTCAAGCGCCTCCGCTCGGTGCCACTAACAAAG
GTGTGTTTTGATTGTGGTGCCAAAAATCCAGCTGGCAAGCATAACCTATGGAGTGTTCCCTTGCATT
GATTGCTCAGGGTCCACCGGTCACTTGGTGTTCACCTGAGTTTTATTCCATCTACAGAGTTGGATTCC
AACTGGTCATGGTTTCAGTTGCGATGCATGCAAGTCGGAGGAAACGCTAGTGCACGTGGCTTGATAGT
TGTGTGGTTCCACCTTTGTCCCCTCCACCAAAGGAGGAAGATTTTTTGCCTCTCACGTTTCTCCTGAG
GTGAGTGACACAGCGTGGGCATCAGCAATAGCAGAACCATCTTCTTAAACATCAAGGCCTGTGAAACC
ACTTTGAAAAATAATGAAGGTGGACAAGAGCAAGGACCAAGTGTGGAAGTCTTAATGTACCAACAAAG
GCTACTTTAGAGGTATCCTCTATCATAAAAAAGAAACCAATCAAGCTAAAAAGGCCTTGGGGCCAAA
AAAGGAAGTTGGGAGCTCAGAACTGGCAAACACATGCTTTAATGAAATTGAAAAACAAGCTCAAGCT
GCGGATAAAATGAAGGAGCAGGAAGACCTGGCCAAGGTGGTATCTAAAGAAGAATCAATTGTTTCATCA
TTACGATTAGCCTATAAGGATCTTGAAATTCAAATGAAGAAAGACGAAAAGATGAACATTAGTGGCAAA
AAAAATGTTGACTCAGACAGACTCGGCATGGGATTTGAAATTCAGAAAGTGTATTTCACATTCAGTG
ACTTCAGATATGCAGACCATAGAGCAGGAATCACCCATTATGGCAAAACCAAGAAAAAGTATAATGAT
GACAGTGACGATTCAATTTTACTTCCAGCTCAAGTTACTTTGACGAGCCAGTGGAGTTAAGGAGCAGT
TCTTTCTCTAGCTGGGATGACAGTTCAGATTCCTATTGGAAAAAGAGACCAGCAAGATACTGAAACA
GTTCTGAAACCACAGGCTATTCAGACAGACCTACTGCTCGCCGAAGCCAGATTATGAGCCAGTTGAA
AATACAGATGAGGCCAGAAGAAGTTGGCAATGTCAAGGCCATTCATCAGATATGATATTTTGAAGA
CAATCCCAGGCTGATTATGAGACCAGGCCCGCCTAGAGAGGCTGTCGGCAAGTTCCTCCATAAGCTCG
GCTGATCTGTTTCGAGGAGCCGAGGAAGCAGCCAGCAGGGAACACAGCCTGTCCAGTGTGCTGCCAAC
GCCCCGACATGGCGCAGTTCAAGCAGGGAGTGAGATCGTTGCTGAAAACTCTCCGTCTTTGCTAAT
GGAGTCGTGACTTCAATTCAGGATCGCTACGGTCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001142293

Insert Size: 1419 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142293.1</u>
RefSeq Size:	2719 bp
RefSeq ORF:	1419 bp
Locus ID:	26286
UniProt ID:	<u>Q9NP61</u>
Cytogenetics:	22q13.2
Protein Pathways:	Endocytosis
MW:	52 kDa
Gene Summary:	<p>The protein encoded by this gene is a GTPase-activating protein (GAP) that associates with the Golgi apparatus and regulates the early secretory pathway of proteins. The encoded protein promotes hydrolysis of ADP-ribosylation factor 1 (ARF1)-bound GTP, which is required for the dissociation of coat proteins from Golgi-derived membranes and vesicles. Dissociation of the coat proteins is a prerequisite for the fusion of these vesicles with target compartments. The activity of this protein is sensitive to phospholipids. Multiple transcript variants encoding different isoforms have been found for this gene. This gene was originally known as ARFGAP1, but that is now the name of a related but different gene. [provided by RefSeq, Nov 2008]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 5' coding region, compared to variant 1. The resulting protein (isoform 2) has the same N- and C-termini but is shorter than isoform 1.</p>