

Product datasheet for SC325988

EFCAB3 (NM_001144933) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	EFCAB3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325988 representing NM_001144933. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGGAGTTGATGCCTCATTTGACATTTAACGGAATGGAAGATTAATGTTAATTCAATAATGGAG
GGCCTGAAGAAGTTTAAACAAAAGGAATGGTGACTCTGCATAAGTTGAAACTGCAATGATATTTAAA
GACAGAGTCACTGGCCACATGGCAGTTTCAGAAATTAACCAAACTTAAGCTGAATCCTCTAACAAAA
GTACCATCTCCCAACAATAAAGGGATAGAGACTTACCAGGATCTCTTCAATGCCAATTACAACACAAA
GAAAAGAAGCTAAGTGCTTCACAAATGGCAGCTTTCCAAGATGCCTACAACCTCTTCTACAAGGACAAA
ACTGGCTGTATTGATTTCCATGGACTGATGTGCACTGTAGCTAAGCTGGGAATGAATCTGACCAAGCAT
GATGTCTATAATGAATTGAAATGTGCTGATTTGATCGAGATGGGAAGGTGAACCTCTCAGACTTTATC
AAGGTTCTAACAGATAAGAATCTCTTCTCAAGGCAGTGTTCCAGAAAAGGAGACCTGTTTAGATTG
GCTGGCAACCCAGGAATCCTATTGTTTGAATCCTATCAAGGCTTCTAGAGACTTCAGCCCTACCCAGA
AAGTCTATAATAGAAATAGTAAGCTATTTCCAAGAAAATTCAGCATACTGGCCAGGAATGTTGTGG
AGTCCCTACACTATGGGCTATGAAAAAGGACACTTAAGCCAGACATATGCACACCTCCAAGCTCAAGC
ATGGCTGCCTTTGCTAATGCTGCCCGATTGCAATAATGAAAGAAAAGGATTTATTTAAATTTCTTGAA
GAGCTCAAGAGATGCAATTCGGTTTCAGATAGCCCATATTCAAAAATACCCATCTTTCCATTGTTCCCT
AATGTGGATGGGGTGGTATGGGAAAGCCATTCAAAGACATGCAGAAATAGAGATGCTAAGAATAAAG
GAGCCTTTGCATTTCTTTGAGGATTATTTTTCCATAAAAGAGACTGGAAAACACAGGCAGCAATATA
AAGTCTATGGACCTGCCTCAGGTTATTCAAATAACATCTTCAACATTGATCAAAATGCTCAAGAAAAAG
CAGACTTGACAGTGGCCGATGCAACTGCTATTAACAACATGTGAAGAGAGCTACTGATACCTATAAT
TTAGGAATTGCCCTTGAACACCGAAAAGAGATGCTAAACCTCTGGCAGAAGATCCGAGGGGATTTGATT
GGGATGGACTCTAGAAATGAGTCCTTTTATGACACCTTTTCTACTATACATGGTCTGGAATGTTTGC
CAAGAATTACTTTCTCCTAAGGACTTAAGGTTATATGATGCCTATGTGAATAGAAATTCCTCCCATAAC
TCCAGATCCTCTTCTCATCAGATACCAAGTGAATGTTACACAGACTCAGGAAGAAAAAGAAAACGGAAA
GGTTTAAAGGGATTTCAACAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001144933
Insert Size:	1473 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001144933.1
RefSeq Size:	1694 bp
RefSeq ORF:	1473 bp
Locus ID:	146779
UniProt ID:	Q8N7B9
Cytogenetics:	17q23.2
MW:	56 kDa