

Product datasheet for **SC326005**

IMPDH1 (NM_001142574) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IMPDH1 (NM_001142574) Human Untagged Clone
Tag:	Tag Free
Symbol:	IMPDH1
Synonyms:	IMPD; IMPD1; IMPDH-I; LCA11; RP10; sWSS2608
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326005 representing NM_001142574.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGACTACCTGATCAGCGGGCACCGGCTACGTGCCCGAGGATGGGCTCACCGCGCAGCAGCTC
TTGGCCAGCGCCGACGGCTCACCTACAACGACTTCTGATTCTCCAGGATTCATAGACTTCATAGCT
GATGAGGTGGACCTGACCTCAGCCCTGACCCGGAAGATCACGCTGAAGACGCCACTGATCTCTCCCCC
ATGGACACTGTGACAGAGGCTGACATGGCCATTGCCATGGCTCTGATGGGAGGTATTGGTTTCATTAC
CACAACCTGCACCCAGAGTTCCAGGCCAACGAGAAGTTTGAACAGGGCTTCATCACGGACCTGTGGTG
CTGAGCCCTCGCACACTGTGGCGATGTGCTGGAGGCCAAGATGCGGCATGGCTTCTCTGGCATCCCC
ATCACTGAGACGGGCACCATGGGCAGCAAGCTGGTGGGCATCGTCACCTCCCAGACATCGACTTCTT
GCTGAGAAGGACCACACCACCTCCTCAGTGAGGTGATGACGCCAAGGATTGAAGTGGTGGTGGCTCCA
GCAGGTGTGACGTTGAAAGAGGCAAATGAGATCCTGCAGCGTAGCAAGAAAGGGAAGCTGCCTATCGTC
AATGATTGCGATGAGCTGGTGGCCATCATCGCCCGCACCGACCTGAAGAAGAACCAGACTACCCTCTG
GCCTCAAAGGATCCCAGAAGCAGCTGCTCTGTGGGGCAGCTGTGGGCACCCGTGAGGATGACAAATAC
CGTCTGGACCTGCTCACCGAGCGGGCGTGCAGCTCATAGTCTTGGACTCGTCCCAGGGAATTCGGTG
TATCAGATCGCCATGGTGCATTACATAAACAAGTACCCCCACCTCCAGGTGATTGGGGGAACGTG
GTGACAGCAGCCCAGGCCAAGAACCTGATTGATGCTGGTGTGGACGGGCTGCGCGTGGGCATGGGCTGC
GGCTCCATCTGCATCACCCAGGAAGTATGGCCTGTGGTGGCCCCAGGGCACTGCTGTGTACAAGGTG
GCTGAGTATGCCCGGCGCTTTGGTGTGCCATCATAGCCGATGGCGGCATCCAGACCGTGGGACACGTG
GTCAAGGCCCTGGCCCTTGGAGCCTCCACAGTGATGATGGGCTCCCTGCTGGCCGCCACTACGGAGGCC
CCTGGCGAGTACTTCTTCTCAGACGGGGTGGGCTCAAGAAGTACCGGGGCATGGGCTCACTGGATGCC
ATGGAGAAGAGCAGCAGCAGCCAGAAACGATACTTCAGCGAGGGGATAAAGTGAAGATCGCGCAGGGT
GTCTCGGGCTCCATCCAGGACAAAGGATCCATTAGAAGTTCGTGCCCTACCTCATAGCAGGCATCCAA
CACGGCTGCCAGGATATCGGGGCCCGCAGCCTGTCTGCTTCCGTCCATGATGTAATCAGGAGAGCTC
AAGTTTGAGAAGCGGACCATGTCGGCCAGATTGAGGGTGGTGCATGGCTGCACTCTTACGAAAAG
CGGCTGTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001142574

Insert Size:	1530 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142574.1
RefSeq Size:	2380 bp
RefSeq ORF:	1530 bp
Locus ID:	3614
UniProt ID:	P20839
Cytogenetics:	7q32.1
Protein Families:	Druggable Genome
Protein Pathways:	Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism
MW:	54.8 kDa
Gene Summary:	<p>The protein encoded by this gene acts as a homotetramer to regulate cell growth. The encoded protein is an enzyme that catalyzes the synthesis of xanthine monophosphate (XMP) from inosine-5'-monophosphate (IMP). This is the rate-limiting step in the de novo synthesis of guanine nucleotides. Defects in this gene are a cause of retinitis pigmentosa type 10 (RP10). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR and coding sequence and lacks an alternate in-frame segment compared to variant 1. The resulting isoform (f) is shorter at the N-terminus and lacks an internal segment compared to isoform a.</p>