

Product datasheet for **SC326006**

DKC1 (NM_001142463) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DKC1 (NM_001142463) Human Untagged Clone
Tag:	Tag Free
Symbol:	DKC1
Synonyms:	CBF5; DKC; DKCX; NAP57; NOLA4; XAP101
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC326006 representing NM_001142463.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGGATGCGGAAGTAATTATTTTCCAAAGAAACATAAGAAGAAAAAGGAGCGGAAGTCATTGCCA
GAAGAAGATGTAGCCGAAATACAACACGCTGAAGAATTTCTTATCAAACCTGAATCCAAAGTTGCTAAG
TTGGACACGCTCTCAGTGGCCCTTTTGCTAAAGAATTTTGATAAGCTGAATGTAAGGACAACACACTAT
ACACCTCTTGCATGTGGTTCAAATCCTCTGAAGAGAGAGATTGGGACTATATCAGGACAGGTTTCATT
AATCTTGACAAGCCCTCTAACCCCTCTCCCATGAGGTGGTAGCCTGGATTTCGACGGATACTTCGGGTG
GAGAAGACAGGGCACAGTGGTACTCTGGATCCCAAGGTGACTGGTTGTTAATCGTGTGCATAGAACGA
GCCACTCGCTTGGTGAAGTCACAACAGAGTGCAGGCAAAGAGTATGTGGGGATTGTCGGCTGCACAAT
GCTATTGAAGGGGGACCCAGCTTTCTAGGGCCCTAGAACTCTGACAGGTGCCTTATCCAGCGACCC
CCACTTATTGCTGCAGTAAAGAGGCAGTCCGAGTGAGGACCATCTACGAGAGCAAAATGATTGAATAC
GATCCTGAAAGAAGATTAGGAATCTTTTGGGTGAGTTGTGAGGCTGGCACCTACATTCGGACATTATGT
GTGCACCTTGTTTGTATTGGGAGTTGGTGGTCAGATGCAGGAGCTTCGGAGGGTTCGTTCTGGAGTC
ATGAGTAAAAGGACCACATGGTGACAATGCATGATGTGCTTGATGCTCAGTGGCTGATGATAACCAC
AAGGATGAGAGTTACCTGCGCGAGTTGTTTACCCTTTGGAAAAGCTGTTGACATCTCATAAACGGCTG
GTTATGAAAGACAGTGCAGTAAATGCCATCTGCTATGGGGCCAAGATTATGCTTCCAGGTGTTCTTCGA
TATGAGGACGGCATTGAGGTCAATCAGGAGATTGTGGTTATCACCACAAAGGAGAAGCAATCTGCATG
GCTATTGCATTAATGACCACAGCGGTCTCTCTACCTGCGACCATGGTATAGTAGCCAAGATCAAGAGA
GTGATCATGGAGAGAGACTTACCCTCGAAGTGGGGTTTAGGTCCAAGGCAAGTCAGAAGAAGCTG
ATGATCAAGCAGGCCCTTCTGGACAAGCATGGGAAGCCACAGACAGCACACCTGCCACCTGGAAGCAG
GATGAGTCTGCCAAAAAGAGGTGGTTGCTGAAGTGGTAAAAGCCCGCAGGTAGTTGCCGAAGCAGCA
AAAAGTGCGAAGCGAAGCGAGAGAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT
ATCAAGAAGGAAAAGAAGAGTAAGAAGGACAAGAAGGCCAAAGCTGGTCTGGAGAGCGGGGCCGAG
CCTGGAGATGGGGACAGTGATACCACCAAGAAGAAGAAGAAGAAGAAGAAGAAGCAAAAGAGGTAGAATTG
GTTTCTGAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001142463

Insert Size:	1530 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142463.2
RefSeq Size:	2593 bp
RefSeq ORF:	1530 bp
Locus ID:	1736
UniProt ID:	O60832
Cytogenetics:	Xq28
Protein Families:	Druggable Genome
MW:	57 kDa
Gene Summary:	<p>This gene functions in two distinct complexes. It plays an active role in telomerase stabilization and maintenance, as well as recognition of snoRNAs containing H/ACA sequences which provides stability during biogenesis and assembly into H/ACA small nucleolar RNA ribonucleoproteins (snoRNPs). This gene is highly conserved and widely expressed, and may play additional roles in nucleo-cytoplasmic shuttling, DNA damage response, and cell adhesion. Mutations have been associated with X-linked dyskeratosis congenita. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. The resulting isoform (2) is shorter than isoform 1.</p>