

## Product datasheet for **SC326010**

### ETV3 (NM\_001145312) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ETV3 (NM_001145312) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETV3
Synonyms:	METS; PE-1; PE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC326010 representing NM\_001145312.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAAAGCCGGCTGTAGCATCGTGGAAAAGCCAGAAGGAGGTGGAGGGTATCAGTTTCTGACTGGGCC
TACAAAACAGAGTCATCCCAGGCTCCCGCAGATCCAGCTGTGGCACTTCATCCTGGAGCTGCTGCAG
AAGGAAGAGTTCCGCCATGTCATCGCTGGCAGCAGGGAGAGTACGGGAATTTGTCATCAAGGATCCA
GATGAGGTGGCCCGCTCTGGGGCCGCAGGAAATGCAAACCACAGATGAATTATGACAAGCTGAGCCGG
GCCCTCAGATACTATTACAACAAGAGGATCCTTCATAAAACAAAAGGAAAAGATTTACCTATAAATTT
AACTTCAACAAGCTGGTATGCCAACTACCCATTTCATCAACATTCGGTCAAGTGGTGTGGTTTCTCAG
AGTGCACCACCAGTGCCAACAGCCTCTTCCCGTTCCATTTCCACCTCTGGACACCCATTCTCCAACC
AATGATGTGCAGCCGGACGGTTCTCTGCTAGCTCCCTAACTGCTTCTGGCCAGGAGTCCAGTAATGGT
ACTGATAGAAAGACTGAGCTTTCAGAGCTGGAGGATGGCTCAGCTGCTGACTGGCCCGGGGTGTGGAT
CCCGTGTCTCCAGGAATGCCATTGGTGGAGGAGGGATTGGCCATCAGAAACGCAAGCCTGACATAATG
CTTCTCTGTTTCTAGGCCAGGGATGTACCCTGACCCCCACAGTCCCTTCGCTGTCTCTCAATCCCT
GGCCGGGAGGTGTCTTAATGTCCCCATTTACCAGCCCTCTCCCTGACTCCCACCATCTTCTCTAT
AGCCCATCACCAGGCCTGAGCCCTTTCACCAGCAGCAGCTGCTTCTCTTCAACCCTGAGGAAATGAAA
CACTACCTTATTCTCAAGCTTGTCTGTGTTCAACTACCATCTGAGTCCTCGGACTTTTCCCGTTAC
CCAGGGCTCATGGTTCACCACTGCAGTGCCAAATGCATCCTGAGGAGTCAACTCAGTTCTCCATCAAG
CTGCAGCCCCACCAGTTGGGCGGAAGAACCGGAGAGGGTTGAGAGCAGCGAGGAGTCAAGCACCAGTC
ACCACGCCACCATTGGCTTCTATCCCCAAGAATCAAGTGGAACCTGCCTCTGAAAAGGATCCTGAG
AGCCTCAGGCAGTCGGCAGGAGAAGGAGGAGCAGCACTCAAGAAGAGGGCACTGTGCCAAGCAGGACC
ATTGAAGAGGAAAAAGGCACCATCTTTGCCGTCTCTGTCACCACCCATCTGGCCCTCTGTGCCATT
AGCACCCCAAGTGGAGAACCCTCTGGAGGTGACTGAAGACAGTGAAGATAGGCCTGGCAAAGGCCAGT
GCACCTGAGAAGAAAGATGCACTGATGCCCCCAAGCTTCGGTTGAAGCGGCGCTGGAATGATGAC
CCTGAAGCCCAGAGCTGAGCAAGAGTGGAAGTTTCTCTGGAATGGGTGAGACCCAGGGCTTGGCA
ACAGCAGCTGCTGATGCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_001145312

<b>Insert Size:</b>	1539 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001145312.1</a>
<b>RefSeq Size:</b>	2002 bp
<b>RefSeq ORF:</b>	1539 bp
<b>Locus ID:</b>	2117
<b>UniProt ID:</b>	<a href="#">P41162</a>
<b>Cytogenetics:</b>	1q23.1
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	57 kDa
<b>Gene Summary:</b>	<p>Transcriptional repressor that contribute to growth arrest during terminal macrophage differentiation by repressing target genes involved in Ras-dependent proliferation. Represses MMP1 promoter activity.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and predicted protein isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>