

Product datasheet for **SC326021**

ZNF286 (ZNF286A) (NM_001130842) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ZNF286 (ZNF286A) (NM_001130842) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ZNF286A |
| Synonyms: | ZNF286 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)


Fully Sequenced ORF: >SC326021 representing NM_001130842.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGACAGATTTGGCTGAAATGCCTGAGAAAGGAGCTCTGTCTTCCCAGGATTCTCCCCATTCCAA
GAGAAGACACAGAAGAGGGGAGAAGTGGCTGCTCTGCGCCTCACGGCCAGATCCCAGGAAACAGTGACA
TTCAAGGATGTGGCCATGGACTTTACACCAGAGGAGTGGGGGAAGCTGGATCCTGCACAAAAGGGATGTG
ATGCTGGAGAACTATAGGAACCTAGTCTCACTTTGGCTTCCAGTTTCCAAACCTGAGAGCTACAACCTG
GAGAATGGAAAAGAACCATTGAAGCTTGAGAGAAAAGCCCCAAAAGCAGCTATTCAGACATGGAGACT
AGACCACAGAGCAAGGATTCAACTTCAGTGCAAGATTTTCCAAAGCAGAATCATGCAAAGTTGCAATA
ATAGACAGACTGACACGGAATAGTGTCTATGACTCTAAGTGAAGCAGCCCTTGAATGTGAAAATTGG
TTAGAGAATCAGCAAGGAAATCAGGAGAGACATTTGAGAGAAATGTTCACTCACATGAATCACTCTCT
GAGGAAACAGACCATAAGCATGATGATACTGGAAAAGCTTCAATCAGAAATCTGTCCTTATCACTGAA
GACAGAGTTCCCAAAGGATCTTATGCCCTCCATACACTTGAAAAAGCTTGAAACAAAAATCAAACCTTA
ATGAAAAAGCAGAGAACTTATAAAGAGAAAAACCTCATAAATGTAATGATTGTGGTGAACCTTCCACC
TACCATTCACTGCTTATTCGACACCAGAGAGTCCATACTGGAGAGAAAACCTATACCTGCAATGAATGT
GGGAAATCTTTTAGCCACAGAGCTAATTTAACTAAACACCAGAGAAGCTATACTAGAATTCTTTTGAA
TGCAGTGAATGCAAGAAAACCTTACAGAAAAGCTCATCCCTTGCAACACATCAGAGAATTCACGTTGGA
GAGAGACCTTATGAATGCAATGAATGTGGGAAAGGTTTTAATCGAAGTACACATCTTGTGCAGCATCAG
TTGATTCATACTGGAGTGAAGCCTTATGAATGTAATGAATGTGATAAAGCTTTTATTTCATTCAGCA
CTCATTAACATCAAAGAAGCTCACTGGAGAGAAAACCTTATAAATGTCAAGAATGTGGGAAAGCCTTT
AGCCATTGCTCATCCCTAACTAAGCATCAGAGAGTTCACTACTGGAGAAAAGCCATATGAATGCAGTGAA
TGTGGAAAAAAGCTTTTAGTCAGAGCACACATCTTGTTC AACATCAGAGAATTCACACTGGAGAGAAACCC
TATGAGTGAATGAATGTGGGAAAACCTTCAGCCGGAGCTCCAATTTTGCTAAACATCAAAGAATTCAT
ATTGAAAAGAAACCGTACAAATGTAGCGAGTGTGGAAAAGCCTTCATTTCATTCAGCTCTATTCAA
CATCAGAGAAGCTCATACCGGAGAGAAAACCTTTAGATGTAATGAGTGTGGGAAAAGCTTTAAGTGCAGT
TCATCTCTCATCAGACATCAAAGAGTTCACACTGAAGAGCAACCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAAAGCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001130842

| | |
|-------------------------------|---|
| Insert Size: | 1566 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001130842.1 |
| RefSeq Size: | 5472 bp |
| RefSeq ORF: | 1566 bp |
| Locus ID: | 57335 |
| UniProt ID: | Q9HBT8 |
| Cytogenetics: | 17p12 |
| Protein Families: | Transcription Factors |
| MW: | 60.2 kDa |
| Gene Summary: | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream start codon, compared to variant 3. Variants 1 and 2 encode the same isoform (2), which has a shorter N-terminus, compared to isoform 1. |