

Product datasheet for **SC326024**

AGFG1 (NM_001135189) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGFG1 (NM_001135189) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGFG1
Synonyms:	HRB; RAB; RIP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC326024 representing NM_001135189.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGCCAGCGCGAAGCGGAAGCAGGAGGAGAAGCACCTGAAGATGCTGCGGGACATGACCGGCCTC
CCGACAACCGAAAGTGCTTCGACTGCGACCAGCGCGGCCACCTACGTTAACATGACGGTCGGCTCC
TTCGTGTGTACCTCTGCTCCGGCAGCCTGCGAGGATTAATCCACCACACAGGGTGAAATCTATCTCC
ATGACAACATTCACACAACAGGAAATTGAATTCTTACAAAAACATGGAATGAAGTCTGTAACAGATT
TGGCTAGGATTATTTGATGATAGATCTTCAGCAATTCAGACTTCAGGGATCCACAAAAAGTGAAAGAG
TTTCTACAAGAAAAGTATGAAAAGAAAAGATGGTATGTCCCGCCAGAACAAGCCAAAGTCGTGGCATCA
GTTTCATGCATCTATTTAGGGTCTCTGCCAGTAGCACAAGCAGCACACCTGAGGTCAAACCACTGAAA
TCTCTTTTAGGGGATTCTGCACCAACACTGCACCTAAATAAGGGCACACCTAGTCAGTCCCAGTTGTA
GGTCGTTCTCAAGGCAGCAGCAGGAGAAGAAGCAATTTGACCTTTTAAGTGATCTGGCTCAGACATC
TTTGCTGCTCCAGCTCCTCAGTCAACAGCTACAGCCAATTTTGCTAACTTTGCACATTTCAACAGTCAT
GCAGGTGGAAGTGTGCATCAGTAAATGCTAATTTTGCTCAATTTGATAACTTCCCAAATCCTCCAGT
GCTGATTTTGAACCTTCAATACTTCCCAGAGTCATCAAACAGCATCAGCTGTTAGTAAAGTTTCAACG
AACAAAGCTGTTTACAGACTGCAGACAAATATGCAGCACTTGCTAATTTAGACAATATCTTCAGTGCC
GGGCAAGGTGGTATCAGGGAAGTGGCTTTGGGACCACAGGTAAGCTCCTGTTGGTTCTGTGGTTTCA
GTTCCAGTCAGTCAAGTGCATCTTCAGACAAGTATGCAGCTCTGGCAGAACTAGACAGCGTTTTAGT
TCTGCAGCCACCTCCAGTAAATGCGTATACTTCCACAAGTAAATGCTAGCAGCAATGTTTTGGAACAGTG
CCAGTGGTTGCTTCTGCACAGACACAGCCTGCTTCATCAAGTGTGCCTGCTCCATTTGGAGCTAGCCT
TCCACAAATCCATTTGTTGCTGCTGCTGGTCTTCTGTGGCATCTTCTACAAACCCATTTAGACCAAT
GCCAGAGGAGCAACAGCGCAACCTTTGGCACTGCATCCATGAGCATGCCACAGGATTCGGCACTCCT
GCTCCCTACAGTCTTCCACCAGCTTTAGTGGCAGCTTTCAGCAGCCTGCCTTTCAGCCAAAGCAGCT
TTCCCTCAACAGACAGCTTTTTCTCAACAGCCCAATGGTGCAGGTTTTCAGCATTTGGACAAAACAAG
CCAGTAGTAAACCCTTTTGGTCAAGTTCAGCTGCTGGAGTATCTAGTAATCCTTTTATGACTGTGCA
CCAACAGGACAATTTCCAACAGGAAGCTCATCAACCAATCCTTTCTTAGACGCGT
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001135189

Insert Size:	1569 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135189.1
RefSeq Size:	8564 bp
RefSeq ORF:	1569 bp
Locus ID:	3267
UniProt ID:	P52594
Cytogenetics:	2q36.3
MW:	54.2 kDa
Gene Summary:	<p>The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008]</p> <p>Transcript Variant: This variant (4) lacks two alternate in-frame exons and uses an alternate in-frame splice site in the 3' coding region, compared to variant 1, resulting in a shorter protein (isoform 4), compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.</p>