

Product datasheet for **SC326032**

FSD1L (NM_001145313) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSD1L (NM_001145313) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSD1L
Synonyms:	CCDC10; CSDUFD1; FSD1CL; FSD1NL; MIR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326032 representing NM_001145313.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACTCCAGAAATACTGCTTTAAGGAGAATGAAAACGTTACAGTTGATAAAGCCTGTTTCTGATC
TCTAACATCACTATTGGACCAGAGTCTATTAACCTGCAGCAGGAAGCTCTACAGAGGATCATTTCAACT
CTGGCAAATAAAAATGATGAAATTCAGAACTTTATTGATACACTACATCATACTAAAAGGAGTTCAG
GAAAATTCGTCCAACATACTCTCAGAGTTAGATGAAGAATTTGATAGTTTATACTCTATACTGGATGAA
GTAAAAGAAAGTATGATTAACGTATCAAGCAGGAACAAGCTCGTAAATCCCAAGAGTTACAGAGTCAG
ATTAGTCAATGTAATAATGCCCTGGAGAACTCTGAAGAATTAGAAATTTGCAACAAGGTCATTAGAT
ATAAAGGAACCTGAAGAATTTTCAAAGGCTGCCAGACAGATCAAGGATAGAGTCACAATGGCTTCAGCC
TTTCGCCCTTTCTTGAACCAAAGGTCAGTGACAACATGACTCATTTAATGGTGGATTTCTCACAGGAA
AGACAGATGCTGCAAACCTTGAAGTTTTGCCAGTCCCCAAAGCTCCAGAGATAGATCCAGTAGAGTGT
TTGGTGGCAGATAACTCTGTAACAGTGGCTTGGAGAATGCCAGAAGAAGATAATAAGATTGACCATTTT
ATACTGGAACATAGGAAGACTAATTTTGTAGTCCACGTGTAAGGATGAGCGATGCTGGGAGATA
ATTGATAATATTAAGGGTACTGAATATACACTATCAGGCTTAAAAATTTGATTCAAAGTATATGAATTC
AGAGTGGCAGCTTGTAAACAAGGCTGTGGCTGGAGAGTATTCTGATCCAGTACTCTAGAGACCAAAGCA
CTTAACCTCAATTTGGATAACTCCTCATCCCATTTGAACCTGAAAGTTGAAGATACATGTGTAGAGTGG
GATCCTACTGGAGGAAAAGGTCAGAAAGTAAAATTAAGGAAAAGAGAACAAGGCCAGAAAGTGGTACA
CCATCCCCAAAACGAACATCTGTAGGCTCCAGGCCACCAGCAGTAAGAGGCAGTAGAGATCGTTTTACT
GGAGAATCATACACAGTGTGGGAGACTGCTATTGAAAGTGGACAACATTATTGGGAGGTCAAGGCC
CAGAAGGATTGTAATCCTACAGTGTGGGAGTAGCATAAAAACGTTGGGAAAATTTGACCAATTTGGGA
AAGACAAACACTAGTTGGTGTATCCATGTCAACAACCTGGCTACAAAACACATTTGCAGCAAAGCATAAT
AATAAAGTCAAAGCTTTGGATGTTACTGTTCTGAAAAAATAGGTGATTTTTGTATTTTGTAGGGGT
CAACTTTCATTCTATGATGCAAATCTAAACAGTTGCTATATTCTTTAAGACAAAATTTACTCAGCCA
GTACTACCTGGTTTCATGGTATGGTGTGGTGGACTTTCTTTGAGTACTGGGATGCAGGTTCCAAGTGT
GTGAGAACACTTCAGAAAAGTAAAATGGAATGACTGGTTCAGCTAGCAGCCTGAACAATGTTGTTACT
CAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001145313

Insert Size: 1593 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145313.2](#)

RefSeq Size: 7750 bp

RefSeq ORF: 1593 bp

Locus ID: 83856

UniProt ID: [Q9BXM9](#)

Cytogenetics: 9q31.2

MW: 59.6 kDa