

Product datasheet for **SC326033**

LRRC27 (NM_001143757) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC27 (NM_001143757) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC27
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001143757, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGGGAAGCAGCTCCTACGAAGTTCCTCTGTGGCTGCTGCTGATCTGGAGGAGGGT
GCTGGTCAGACTAGGAGCTTGCCCTGCCACCCCTCCAAAGATGTTCAAGGGTGTGGA
GGCATCATCTTTCTCCTCACCATTGTTAGACTTGAGTGAAAGTGGTCTGTGCCGTTG
GAGGAGGTCTTTAGAATCCCCAGCCTTCAACAATTGCATCTGCAAAGGAATGCCCTGTG
GTGATTCCTCAAGATTTCTTTCAGTTGCTTCCGAACCTGACTTGGCTGGACCTCCGGTAC
AATAGAATTAAGCGCTTCCTTCTGGGATTGGAGCTCACCAGCATTGAAACTTTGCTT
TTAGAAAGAAATCCTATCAAAATGTTACCTGTGGAGCTGGGAGCGTAACACGCTGAAA
GCACTGAACCTAAGACACTGCCCTCTGGAATTCCCTCCTCAGCTCGTTGTGCAGAAGGA
TTGGTGGCTATCCAGCGCTTCCTGCGGATGTGGCAGTAGAACAACCTCTCTCCCCAGAAAT
CCAACCTCTCAAGAGGCTCCACCGTTAGAGAGATGACCTCCGTGACCTCCCGAGCCCA
GGACTGGAGTTGTCTGGAGACCACGCTCTAACCAAGGAGCTGTGAACGCTCAGGACCCA
GAGGGGGCTGTGATGAAAGAGAAGGCCAGCTTTCTCCACCTGTGGAAAAGCCAGACCTG
AGTGAACCTCAGGAAGTCTGCGGACTCCTCAGAGAACTGGCCAGCGAGGAGGAGATCAGG
CGCTTTTGAAGCTGAGGCAGGAGATTGTTGAGCACGTGAAGGCAGACGTTCTGGGAGAT
CAGCTCTTGACGAGGGAATTACCTCCAAATCTCAAGGCGCCTTGAACATTGAGAAAGAA
CTACCAAAGCCAAGACACGTTTTTCAGAAGGAAGACAGCCTCCTCCAGGAGCATCTTACCC
GACCTCTGTGCACCGTACCAAAATGGCGATCCGAGCAAAAAGACTGGAAGAGAGCCGAGCG
GCGGCGCTCCGAGAGCTCCAGGAGAAGCAGGCTCTGATGGAGCAGCAGAGACGAGAGAAA
AGGGCACTGCAGGAGTGGAGAGAGCGAGCCAGAGGATGAGGAAGAGGAAGGAAGAGCTC
AGCAAACCTCTGCCTCCGCGGAGGAGCATGGTGGCATCAAAGATTCCTCTGCCACAGAT
CTGATAGATAACAGGAAAGTACCACTGAATCCGCTGGAAAAATGAAACCAAGCAAGAG
AAATCGCCACAAGCAAGTAAAGAAATGAGTGCCCTGCAGGAGAGAAATTTAGAAGAGAAG
ATAAAACAGCACGTCTCCAAATGCGTGAGCAAGAAGATTCCATGGCCAGGCCCACTG
GAGGAGATGAGGAAGGCTGCCGAGGATCTGGAATTTGCCACAGAGCTACAGGATGAAGTA
TTGAAGCTAAATTTGGGATTAACCTTGAACAAAGATCGTCGACGGCGGCCCTCACTGGA
AACCTTTGCTTGGCCTGCCGGCAGCACAGCCTCAAAATACATTTTTTAACAAAAATAT
GGAGAATCAGGAAATGTTTCGAGATACCAG

```

Restriction Sites:	Please inquire
ACCN:	NM_001143757



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001143757.1, NP_001137229.1</u>
RefSeq Size:	8205 bp
RefSeq ORF:	1593 bp
Locus ID:	80313
UniProt ID:	<u>Q9C0I9</u>
Cytogenetics:	10q26.3