

Product datasheet for **SC326040**

ZNF397 (NM_001135178) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF397 (NM_001135178) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF397
Synonyms:	ZNF47; ZSCAN15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326040 representing NM_001135178.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGTGGAATCTGGAGTGATTTCAACCTGATACCTCAGGATCTCCGGAACAAGAATAACTACTA
GTGAAAGTAGAAGATAACTTTTCTGGGATGAGAAATTTAAGCAGAATGGGAGTACTCAATCCTGCCAA
GAATTGTTTCGTAGCAATTCAGAAAATTTGCTACCAGGAGACACCTGGGCCCGGAGGCTCTGAGC
CGACTCCAGGAACTTTGCTATCAGTGGCTAATGCCAGAGTTGCACACAAAGGAGCAGATCTTAGAACTG
CTGGTACTGGAGCAGTTCTGAGCATTCTGCCTGAGGAGCTGCAGATCTGGGTTGAGCAACATAATCCA
GAAAGCGCGGAGGAAGCTGTGACCCTGTTGGAGGATTTAGAGAGGGAGTTTGATGACCCAGGGCAGCAG
GTCCCAGCTAGTCCACAGGGACCAGCAGTGCCATGGAAGGATTTAACATGTCTCAGAGCATCCCAAGAG
TCAACAGACATCCACCTCCAGCCCTTAAGACACAGCTGAAATCTGGAAACCATGCCTTTCCCTAAA
AGTGATTGTGAGAACAGTGAACAGCAACAAAAGAGGGCATCTCAGAAGAAAAATCACAGGGACTCCCT
CAGGAACCTTCATTTGAGGAATTAGTGAGCATGAAAGCAATTTAGTGTGAAGCAAGGAAGTGCTACA
GGGAGAAAATAAGATCTCCTTCCCAAGGGGCGAGTTTTAGTCAAGTGATCTTCACAAACAAATCTCTA
GGAAAGAGAGACCTTTATGATGAGGCTGAAAGATGCTTGATTCTAACTACAGACTCTATAATGTGTGAG
AAAGTTCTCCAGAAGAGAGACCTTATAGATGTGATGTATGTTGGGCACAGCTTCAAGCAGCATTCTCT
CTAACACAACATCAGAGAATCCATACTGGAGAAAAGCCCTATAAATGTAACCAAGTGTGGGAAGGCCCTT
AGTTTGAGGTCCTATCTTATTATTCATCAGAGAATTCATAGTGGTGAGAAAAGCATATGAATGTAGTGAA
TGTGGGAAAGCTTTCAATCAGAGCTCAGCCCTCATTAGACATCGGAAAATCCATACTGGTGAGAAAAGCT
TGTAATGTAATGAGTGTGGCAAAGCATTGAGTCAAAGTTCATATCTCATTATACATCAAAGAATTCAC
ACTGGTGAGAAAACCTTATGAGTGAATGAATGTGGGAAAACCTTTAGCCAGAGCTCAAACCTCATTAGA
CATCAGCGAATTCACACAGGAGAGAGACCTTATGAATGTAATGAATGTGGAAAAGCTTTGAGGCAGAGC
TCAGAGCTGATTACTCATCAGAGAATACATAGTGGAGAGAAAACCTATGAATGTAGTGAATGTGGAAAA
GCTTTCAGTTTGAGCTCAAACCTTATCAGACATCAGAGAATTCATAGTGGGGAGGAACCTTATCAGTGT
AATGAATGTGGCAAAACCTTTCAAAGGAGCTCAGCCCTTGTTCAGCATCAGAGAATTCATTCTGGGGAT
GAAGCTTATATGTAATGAATGTGGGAAGCTTTCAGGCACAGATCGGTCTTATGCGCCATCAAAGA
GTCCACACTATAAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_001135178
Insert Size:	1605 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135178.2
RefSeq Size:	2644 bp
RefSeq ORF:	1605 bp
Locus ID:	84307
UniProt ID:	Q8NF99
Cytogenetics:	18q12.2
Protein Families:	Transcription Factors
MW:	61.1 kDa
Gene Summary:	<p>This gene encodes a protein with a N-terminal SCAN domain, and the longer isoform contains nine C2H2-type zinc finger repeats in the C-terminal domain. The protein localizes to centromeres during interphase and early prophase, and different isoforms can repress or activate transcription in transfection studies. Multiple transcript variants encoding different isoforms have been found for this gene. Additional variants have been described, but their biological validity has not been determined. [provided by RefSeq, Oct 2009]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1), also known as ZNF397-fu.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>