

Product datasheet for **SC326051**

POTEC (NM_001137671) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POTEC (NM_001137671) Human Untagged Clone
Tag:	Tag Free
Symbol:	POTEC
Synonyms:	A26B2; CT104.6; POTE-18; POTE18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326051 representing NM_001137671.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGTGACTGAGGTTTGTCAATGCCCGCTGCCCTGCTGTGAAGAAGCCATTTCGATCTCAGGAGCAAG
ATGGGCAAGTGGTTTACCACCGCTTCCCCTGCTGCAAGGGGAGCGGCAAGAGCAACATGGGCACTTCT
GGAGACCACGACGACTCCTTTATGAAGATGCTCAGGAGCAAGATGGGCAAGTGGTCCACCCTGCTTCC
CCCTGCTGCAGGGGAGCGGCACGAGCAACGTGGGCACTTCTGGAGACCATGACAACCTTTATGAAG
ACGCTCAGGAGCAAGATGGGCAAGTGGTGTCTACTGCTTCCCCTGCTGCAGGGGAGCGGCAAGAGC
AACGTGGGCGCTTGGGAGACTACGACGACAGCGCCTTCATGGAGCCGAGGTACCACGTCGCTCGAGAA
GATCTGGACAAGCTCCACAGAGCTGCCTGGTGGGTAAGGTCCCCAGAAAGGATCTCATCGTCATGCTC
AGGGACACGGACATGAACAAGAGGGACAAGCAAAAGAGGACTGCTCTACATTTGGCCTCTGCCAATGGA
AATTCAGAAGTAGTAACTCCTGCTGGACAGACGATGTCAACTTAATGTCTTGACAACAAAAAAGG
ACAGCTCTGATAAAGCCGTACAATGCCAGGAAGATGAATGTGTGTTAATGTTGCTGGAACATGGCGCT
GATCAAAATATCCAGATGAGTATGGAAATACCACTCTACACTATGCTGTCCACAATGAAGATAAATTA
ATGGCCAAAGCACTGCTCTTATATGGTGTGATATTGAATCAAAAAACAAGTGTGGCCTCACACCACTT
TTGCTTGGCGTACATGAACAAAAACAGCAAGTGGTGAATTTTTAATCAAGAAAAAGCTAATTTAAT
GCCTTGTAGATATGGAAGAACTGCCCTCATACTTGCTGTATGTTGTGGATCAGCAAGTATAGTCAAT
CTTCTCCTTGAGCAAAATGTTGATGTATCTTCTCAAGTCTATCTGGACAGACGGCCAGAGAGTATGCT
GTTTCTAGTCATCATGTAAATTTGTGAATTACTTTCTGACTATAAAGAAAAACAGATGCTAAAAATC
TCTTCTGAAAACAGCAATCCAGAACAAGACTTAAAGCTAACATCAGAGGAAGAGTCACAAGGCTTAAA
GTCAGTAAAAATAGCCAGCCAGAGAAAAATGTCTCAAGAACCAGAAATAAATAAGGACTGTGATAGAGAG
GTTGAAGAAGAAAATAAAGAAGCATGGAAGTAATCCTGTGGGATTACCAGAAAACCTGACTAATGGTGCC
AGTGTGGCAATGGTGTGATGGATTAATTCACAAAGGAGGAGCAGAAAACCTGAAAATCAGCAATTT
CCTGACTGAGAATGAAGAGTATCACAGTGACGAACAAAATGATACCCGAAACAACTTTCTGAAGAA
CAGAACACTGGAATATCACAAGATGAGATTCTGACTAATAAACAAAAGCAGATAGAAGTGGCTGAAAAG
AAAATGAATTCTGAGCTTCTCTTAGTCATAAGAAAGAAGAAGATCTCTTGCCTGAAAACAGCATGTTG
CAGGAAGAAATTGCCATGCTAATTTCTGGAGACTGGAACTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001137671

Insert Size: 1629 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001137671.1](#)

RefSeq Size: 2243 bp

RefSeq ORF: 1629 bp

Locus ID: 388468

UniProt ID: [B2RU33](#)

Cytogenetics: 18p11.21

MW: 61.2 kDa