

## Product datasheet for **SC326052**

### **ABL2 (NM\_001136001) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ABL2 (NM_001136001) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABL2
Synonyms:	ABLL; ARG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC326052 representing NM\_001136001.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGGCAGCAGGTGGGCCGCTCGGGGAAGCTCCGGGGTCCAGCAGCCTCAGCCCCGGGGATCCGG
GGCAGCAGTGCAGCCAGGCCCTCCGGCCGAGGCGGACCCGGCGGGGCGCACCACAGAGACCGGCTTC
AATATCTTCACCCAGCATGAAGCTTGCATCGTCCCTATGGTTGTGATGTTGAACCCAGGCACTAAAT
GAGGCTATCAGGTGGAGCTCCAAGGAGAACTTGCTCGGAGCCACTGAGAGTGACCCTAATCTCTCGTT
GCACTTTATGATTTGTAGCAAGTGGTATAACACACTCAGCATCACTAAAGGTGAAAAGCTACGAGTC
CTTGGTTACAACCAGAATGGTGAGTGGAGTGAAGTTTCGCTCTAAGAATGGGCAGGGCTGGTGCCAAGC
AACTACATCACCCAGTGAACAGCCTGGAAAAACTCCTGGTACCATGGACCTGTGTACGCGAGTGCA
GCTGAGTATCTGCTCAGCAGTCTAATCAATGGCAGCTTCTGGTGCAGAAAAGTGAAGTAGCCCTGGG
CAGCTGTCCATCTCGCTCAGGTACGAGGACGTGTGTACTACTACAGGATCAATACCACTGCAGATGGC
AAGGTGTATGTACTGCTGAGAGCCGCTTCAGCACCTTGGCAGAGCTTGTACACCATCACTCCACAGTG
GCTGATGGGCTGGTGACAACATTACACTACCCAGCACCAAGTGTAAAGCCTACAGTCTATGGTGTG
TCCCCATCCACGACAAATGGGAAATGGAGCGAACAGATATTACCATGAAGCACAAACTTGGGGCGGT
CAGTATGGAGAGGTTTACGTTGGCGTCTGGAAGAAATACAGCCTTACAGTTGCTGTGAAAACATTGAAG
GAAGATACCATGGAGGTAGAAGAATTCCTGAAAGAAGCTGCAGTAATGAAGGAAATCAAGCATCCTAAT
CTGGTACAACCTTTAGGTGTGTACTTTGGAGCCACCATTTTACATTGTGACTGAATACATGCCATAC
GGGAATTTGCTGGATTACCTCCGAGAATGCAACCGAGAAGAGGTGACTGCAGTTGTGCTGCTACATG
GCCACTCAGATTTCTTCTGCAATGGAGTACTTAGAGAAGAAGAATTTTATCCATAGAGATCTTGCAGCT
CGTAAGTGCCTAGTGGGAGAAAACCATGTGGTAAAAGTGGCTGACTTTGGCTTAAGTAGATTGATGACT
GGAGACACTTATACTGCTCATGCTGGAGCCAAATTTCTATTAAGTGGACAGCACCCAGAGAGTCTTGCC
TACAATACCTTCTCAATTAATCTGACGCTGGGCTTTTGGGGTATTGTTGTGGGAAATTGCTACCTAT
GGAATGTCACCATATCCAGGTATTGACCTGTCTCAGGTCTATGACCTACTAGAAAAGGATATCGAATG
GAACAGCCTGAGGGATGCCCCCTAAGGTTTATGAACTTATGAGAGCATGCTGGAAGTGGAGCCCTGCC
GATAGGCCCTCTTTGCTGAAACACACCAAGCTTTTGAACCATGTTCCATGACTCCAGCATTCTGAA
GTTCTCTACATTGTGCCAATCAGACGTGCATAACATTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001136001

**Insert Size:** 1629 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001136001.1</a>
<b>RefSeq Size:</b>	2140 bp
<b>RefSeq ORF:</b>	1629 bp
<b>Locus ID:</b>	27
<b>UniProt ID:</b>	<a href="#">P42684</a>
<b>Cytogenetics:</b>	1q25.2
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	ErbB signaling pathway, Viral myocarditis
<b>MW:</b>	60.9 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the Abelson family of nonreceptor tyrosine protein kinases. The protein is highly similar to the c-abl oncogene 1 protein, including the tyrosine kinase, SH2 and SH3 domains, and it plays a role in cytoskeletal rearrangements through its C-terminal F-actin- and microtubule-binding sequences. This gene is expressed in both normal and tumor cells, and is involved in translocation with the ets variant 6 gene in leukemia. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Nov 2009]</p> <p>Transcript Variant: This variant (e) lacks an alternate in-frame exon in the 5' coding region, and differs in the 3' coding region and 3' UTR, compared to variant b. The resulting isoform (e) has a distinct C-terminus and is shorter than isoform b.</p>