

Product datasheet for **SC326060**

MYH (MUTYH) (NM_001128425) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MYH (MUTYH) (NM_001128425) Human Untagged Clone
Tag:	Tag Free
Symbol:	MUTYH
Synonyms:	MYH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326060 representing NM_001128425.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGACACCGCTCGTCTCCCGCCTGAGTCGTCTGTGGGCCATCATGAGGAAGCCACGAGCAGCCGTGGGA
AGTGGTACACAGGAAGCAGGCAGCCAGCCAGGAAGGGAGGCAGAAGCATGCTAAGAACAACAGTCAGGCC
AAGCCTTCTGCCTGTGATGCCTGTGCAGGGATGATTGCTGAGTGTCTGGGGCCCCAGCAGGCCTGGCC
AGGCAGCCGGAAGAGGTGGTATTGCAGGCCTCTGTCTCCTCATACCATCTATTACAGAGCGTAGCTGAA
GTCACAGCCTTCCGAGGGAGCCTGCTAAGCTGGTACGACCAAGAGAAAACGGGACCTACCTGGAGAAGA
CGGGCAGAAGATGAGATGGACCTGGACAGGCGGCATATGCTGTGTGGGTCTCAGAGGTATGCTGCAG
CAGACCCAGGTTGCCACTGTGATCAACTACTATACCGGATGGATGCAGAAGTGGCCTACACTGCAGGAC
CTGGCCAGTGCTCCCTGGAGGAGGTGAATCAACTCTGGGCTGGCCTGGGCTACTATTCTCGTGGCCGG
CGGCTGCAGGAGGAGCTCGGAAGGTGGTAGAGGAGCTAGGGGGCCACATGCCACGTACAGCAGAGACC
CTGCAGCAGCTCCTGCCTGGCGTGGGGCGCTACACAGCTGGGGCCATTGCCTCTATCGCCTTTGGCCAG
GCAACCGGTGTGGTGGATGGCAACGTAGCACGGGTGCTGTGCCGTGCCGAGCCATTGGTGTGATCCC
AGCAGCACCCTTGTTCCAGCAGCTCTGGGGTCTAGCCCAGCAGCTGGTGGACCCAGCCCGGCCAGGA
GATTTCAACCAAGCAGCCATGGAGCTAGGGGCCACAGTGTGTACCCACAGCGCCCACTGTGCAGCCAG
TGCCCTGTGGAGAGCCTGTGCCGGGCACGCCAGAGAGTGGAGCAGGAACAGCTCTTAGCCTCAGGGAGC
CTGTGGGCAGTCTGACGTGGAGGAGTGTGCTCCCAACTGGACAGTGCCACCTGTGCCTGCCTCCC
TCGGAGCCCTGGGACCAGACCTGGGAGTGGTCAACTTCCCAGAAAGGCCAGCCGCAAGCCCCCAGG
GAGGCCAAGCTCTGCCACCTGTGTTCTGGAACAGCCTGGGGCCCTTGGGGCCAAATCTGCTGGTGCAG
AGTCCCAACTCAGGTCTGCTGGCAGGACTGTGGGAGTCCCGTCCGTGACCTGGGAGCCCTCAGAGCAG
CTTCAGCGCAAGGCCCTGCTGCAGGAACACAGGTTGGGCTGGGCCCTCCAGCCACGCACCTCCGG
CACCTTGGGGAGGTTGTCCACACCTTCTCTCACATCAAGCTGACATATCAAGTATATGGGCTGGCCTTG
GAAGGGCAGACCCAGTGACCACCGTACCACAGGTGCTCGTGGCTGACGCAGGAGGAATTTACACC
GCAGCTGTTTCCACGCCATGAAAAAGTTTTCCGTGTGTATCAGGGCCACAGCCAGGGACCTGTATG
GGTTCCAAAAGTCCCAGGTGCTCTCCGTGCAGTCGAAAAAGCCCCGCATGGGCCAGCAAGTCTCG
GATAATTTCTTTCGGTCTCACATCTCCACTGATGCACACAGCCTCAACAGTGCAGCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001128425
Insert Size:	1650 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001128425.1</u>
RefSeq Size:	1945 bp
RefSeq ORF:	1650 bp
Locus ID:	4595
UniProt ID:	<u>Q9UIF7</u>
Cytogenetics:	1p34.1
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Base excision repair
MW:	60.3 kDa
Gene Summary:	This gene encodes a DNA glycosylase involved in oxidative DNA damage repair. The enzyme excises adenine bases from the DNA backbone at sites where adenine is inappropriately paired with guanine, cytosine, or 8-oxo-7,8-dihydroguanine, a major oxidatively damaged DNA lesion. The protein is localized to the nucleus and mitochondria. This gene product is thought to play a role in signaling apoptosis by the introduction of single-strand breaks following oxidative damage. Mutations in this gene result in heritable predisposition to colorectal cancer, termed MUTYH-associated polyposis (MAP). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2017] Transcript Variant: This variant (alpha5) represents the longest transcript and encodes the longest isoform (5).