

Product datasheet for **SC326064**

RCOR3 (NM_001136223) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RCOR3 (NM_001136223) Human Untagged Clone
Tag:	Tag Free
Symbol:	RCOR3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326064 representing NM_001136223.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCCGGCATGATGGAGAAAGGGCCCGAGTTACTGGGAAGAACCGATCGGCCAACGGCAGCGCCAAG
AGCCCGGCAGCGCGCGCGCAGCGCGCCTCGTCCACCAACGGCGGGCTGCACTACTCAGAGCCCGAG
AGCGGCTGCAGCAGCGACGACGACGACGATGTTGGGATGAGAGTCGGAGCCGAATACCAAGCTCGGATC
CCTGAATTTGATCCAGGTGCTACAAAGTACACAGATAAAGACAATGGAGGGATGCTTGTATGGTCTCCA
TATCACAGTATCCAGATGCCAAATTGGATGAATACATTGCAATTGCAAAGGAAAAGCATGGCTACAAT
GTGGAACAGGCACCTTGGCATGTTGTTCTGGCATAAACATAACATTGAGAAGTCCCTTGCTGATCTCCCT
AATTTCACTCCCTTCCGGATGAGTGGACAGTGGAAAGATAAAGTCCTATTTGAACAAGCCTTTAGTTTT
CATGGAAGAGCTTTACAGGATTCAGCAAATGCTTCCAGATAAGACAATTGCAAGCCTTGTAATAAT
TACTATTCTTGGAAAAAACTCGCTCTAGGACAAGTTTGATGGATCGCCAGGCTCGTAACTAGCTAAT
AGACATAATCAGGGTGACAGTGTATGATGTAGAGAACAACATCCAATGGATGGGAATGATAGTGAT
TATGATCCAAAAAAGAAGCCAAAAAAGAGGGTAATACTGAACAACCTGTCCAACTAGCAAGATTGGA
CTTGGAAGAAGAGAGTATCAGAGTTTACAACATCGCCATCATTCTCAGCGTTCTAAGTGCCGTCCACCT
AAGGGCATGTATTTAACCAGGAAGATGTGGTAGCAGTTTCTGTAGTCCCAATGCGACCCAACACCATC
CTGAGGCAACTGGACATGGAGTTGATCTCTCTAAAACGTCAGGTTTCCAGAAATGCTAAGCAAGTAAACAGT
GCATTTAAACAGAAAAATGGAAGTTGGAATTGAAGAATTCAAACCTCCTGAGTCAAATCAGAAAAATTAAT
GCCCGTTGGACCACAGAGGAGCAGCTTCTAGCAGTGAAGGTGTCGCAAAATATGGTAAAGATTTCAA
GCTATTGCAGATGTAATTGGCAACAAGACTGTTGGCCAAGTGAAGAACTCTTTGTAACACAGGCGT
CGGTTTAACTTAGAGGAGGTATTGCAGGAGTGGGAAGCAGAACAAGGAACCCAGGCTTCTAATGGTGAT
GCTTCTACTTTTAGGGGAGGAGACAAAAAGTGCTTCTAATGTGCCATCAGGGAAGAGCACTGATGAAGAA
GAGGAGGCACAGACCCACAGGCTCCTCGGACACTGGGTCCATCACCTCCTGCCCATCATCCACTCCA
ACACCAACAGCCCTATTGCCACTCTGAACCAGCCTCCACCATTCTTCGTCCAACACTGCCTGCTGCC
CCGGCTCTTACCAGGAGCCTCCTCCACTCCAGCAGCAGGCTCGGTTTATCCAGCCCCGGCCAACTTTA
AATCAGCCTCCACCACCTCTTATTCGCCCTGCTAATTCATGCCACCCGCTCAAACCAAGACCGGTG
TTGTCCACGGTTGGTGGTCAACAGCCACCATCACTTATTGGAATTCAGACAGATTCACAGTCTCACTG
CACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_001136223
Insert Size:	1662 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001136223.1</u>
RefSeq Size:	4474 bp
RefSeq ORF:	1662 bp
Locus ID:	55758
UniProt ID:	<u>Q9P2K3</u>
Cytogenetics:	1q32.2-q32.3
Protein Families:	Transcription Factors
MW:	61.2 kDa
Gene Summary:	<p>May act as a component of a corepressor complex that represses transcription. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.</p>