

## Product datasheet for **SC326066**

### **RNF156 (MGRN1) (NM\_001142289) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	RNF156 (MGRN1) (NM_001142289) Human Untagged Clone
Tag:	Tag Free
Symbol:	MGRN1
Synonyms:	RNF156
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC326066 representing NM\_001142289.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCTCCATTCTCAGCCGCCGCATCGCGGGGTGGAGGACATCGACATCCAGGCGAACTCGGCCTAT
CGCTACCCTCCGAAGTCCGGAACTACTTTGCTTCGCACTTTTTCATGGGAGGAGAGAAATTCGACACC
CCCCACCCTGAAGGTTACCTCTTTGGAGAGAACATGGATCTGAACCTTCTGGGCAGCCGCCGGTCCAG
TTTCCCTACGTCACTCCTGCCCCACGAGCCGTGAAGACGCTGCGGAGCCTGGTGAACATCCGCAAA
GACTCCCTGCGGCTGGTGAAGTACAAAGACGATGCCGACAGCCCAACGAGGACGCGACAAGCCCCGG
GTGCTCTACAGCCTGGAGTTCACCTTCGACGCCGATGCCCGCGTGGCCATCACCATCTACTGCCAGGCA
TCGGAGGAGTTCTGAACGGCAGGGCAGTATACAGCCCCAAGAGCCCTCGCTACAGTCCGAGACCGTC
CACTACAAGAGAGGGGTGAGCCAGCAGTCTCCCTGCCCTCCTCAAGATTGACTTCTCGGAATGGAAG
GATGACGAGCTGAACTTTGACCTGGACCGGGCGTGTTCAGTAGTCATCCAGGCTGTGGTGGACGAA
GGAGATGTGGTGAAGTACTGGCCACGCCACGTGCTTTGGCTGCCTTTGAAAAGCACATGGACGGC
AGCTTCTCTGTGAAGCCTTTAAAGCAGAAGCAAATTTGGACCGGTGAGTACCTCCTGCAGGAGATC
TATGGCATTGAGAACAAGAACAACCAGGAGACCAAGCCCTCGGACGACGAGAACAGCACAACAGCAAC
GAGTGTGTGGTGTGCCTGTCCGACCTGCCGGACACGCTGATCCTGCCCTGCCGCCACCTGTGCCTGT
ACCTCCTGCGCCGACACGCTGCGCTACCAGGCCAACAAGTCCCATCTGCCGGTGCCTTTCCGGGCC
CTCCTGCAGATCCGGGCGGTGCGGAAGAAGCCAGGAGCCCTGTCCCCGTGTCTTCAGCCCCGTCTG
GCCCAGAGCCTGGAGCATGATGAGCACTAATCTGACAGCGTCCCACCTGGCTACGAGCCCATCTCG
CTGCTCGAGGCGCTCAACGGCCTCCGGGCTGTCTCCCCGGCCATCCCCTCGGCCCTCTTATGAAGAA
ATCACCTATTCAGGCATCTCGGACGGCCTGTCCCAGGCCAGCTGTCCCCTCGCGGCTATCGACCACATC
CTGGACAGCAGCCGCCAGAAGGGCAGGCCGAGAGCAAGGCCCCGACAGCACCTACGGTCCCCGTCT
TCCCCATCCACGAAGAGGATGAGGAGAAGCTCTCCGAGGACGTGGACGCCCTCCCCACTGGGTGGC
GCAGAGCTGGCCCTGCGGAAAGCAGCTCCCCTGAGAGTTTCATAACAGAAGAGGTTGATGAGTCGTG
TCACCACAGCAAGGGACCCGAGCAGTTCATTGAGAATGTCTGCAGGACAGCAGCCCCGAGCACTGT
GGCCGAGGCCACCTGCTGACATCTACCTGCCAGGACGGCCACCTCCATGGAGACGGCCACGGCCTC
GCCACCACGAGCCACCTGGCCTCCACTTGGTGGCCCCAGCCCCGATCCAGCGCCGCGAGCTGACC
CCACTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001142289

**Insert Size:** 1665 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001142289.2](#)

**RefSeq Size:** 3868 bp

**RefSeq ORF:** 1665 bp

**Locus ID:** 23295

**UniProt ID:** [O60291](#)

**Cytogenetics:** 16p13.3

**Protein Pathways:** Ubiquitin mediated proteolysis

**MW:** 60.7 kDa

**Gene Summary:** Mahogunin (MGRN1) is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in vitro.[supplied by OMIM, Apr 2004]  
Transcript Variant: This variant (2) lacks an alternate in-frame exon, compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.  
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.