

Product datasheet for **SC326072**

IMPDH1 (NM_001142576) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IMPDH1 (NM_001142576) Human Untagged Clone
Tag:	Tag Free
Symbol:	IMPDH1
Synonyms:	IMPD; IMPD1; IMPDH-I; LCA11; RP10; sWSS2608
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326072 representing NM_001142576.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGGGCCACTCACTCCACCACCGCTGCAGGGAGGCGGAGCCGCGCTGTTCCGGAGCCCGGAGCC
CGGCAACACCCGGGACACGAGACGGCGGGCAGCGGTACAGCGCCCGACTGCTGCAGGCCGGCTACGAG
CCCGAGAGCCCTAGATTGGACCTCGCTACACACCCGACGACACCCCGTTCAAGAACTATCTTCAGTGGTC
TACTGGCAGGTGTTGGTGTCCAGATGGATCGCCTTCGACAGGGCTAGCGACTTCTGATTCTCCAGGA
TTCATAGACTTCATAGCTGATGAGGTGGACCTGACCTCAGCCCTGACCCGGAAGATCACGCTGAAGACG
CCACTGATCTCTCCCCATGGACTGTGACAGAGGCTGACATGGCCATTGCCATGGCTCTGATGGGA
GGTATTGGTTTCATTACCACAACCTGCACCCAGAGTTCCAGGCCAACGAGGTGCGGAAGGTCAAGAAG
TTTGAACAGGGCTTCATACGGACCCTGTGGTGTGAGCCCTCGCACACTGTGGCGATGTGCTGGAG
GCCAAGATGCGGCATGGCTTCTCTGGCATCCCCATCACTGAGACGGGCACCATGGGCAGCAAGCTGGTG
GGCATCGTCACCTCCCGAGACATCGACTTCTTGTGCTGAGAAGGACCACACCACCCTCCTCAGTGAGGTG
ATGACGCCAAGGATTGAACTGGTGGTGGCTCCAGCAGGTGTGACGTTGAAAGAGGCAAATGAGATCCTG
CAGCGTAGCAAGAAGGGGAAGCTGCCTATCGTCAATGATTGCGATGAGCTGGTGGCCATCATCGCCCGC
ACCGACCTGAAGAAGAACCGAGACTACCTCTGGCCTCCAAGGATTCCCAGAAGCAGCTGCTCTGTGGG
GCAGCTGTGGGCACCCGTGAGGATGACAAATACCGTCTGGACCTGCTCACCCAGGCGGGCGTCGACGTC
ATAGTCTTGACTCGTCCAAGGGAATTCGGTGTATCAGATCGCCATGGTGCATTACATCAAACAGAAG
TACCCACCTCCAGGTATTGGGGGAACGTGGTGACAGCAGCCAGGCCAAGAACCTGATTGATGCT
GGTGTGGACGGCTGCGCGTGGGCATGGCTGCGGCTCCATCTGCATCACCCAGGAAGTATGGCCCTGT
GGTCGGCCCCAGGGCACTGCTGTGTACAAGGTGGCTGAGTATGCCCGGCGCTTTGGTGTGCCCATCATA
GCCGATGGCGGCATCCAGACCGTGGGACACGTGGTCAAGGCCCTGGCCCTTGAGCCTCCACAGTGATG
ATGGGCTCCCTGCTGGCCGCACTACGGAGGCCCTGGCGAGTACTTCTTCTCAGACGGGTGCGGCTC
AAGAAGTACCGGGCATGGGCTCACTGGATGCCATGGAGAAGAGCAGCAGCAGCCAGAAACGATACTTC
AGCGAGGGGGATAAAGTGAAGATCGCGCAGGGTGTCTCGGGCTCCATCCAGGACAAGGATCCATTAG
AAGTTCGTGCCCTACCTCATAGCAGGCATCCAACACGGCTGCCAGGATATCGGGGCCCGCAGCCTGTCT
GTCCTTCGGTCCATGATGTACTCAGGAGAGCTCAAGTTTGAGAAGCGGACCATGTCGGCCAGATTGAG
GGTGGTGTCCATGGCCTGCACTTTACGAAAAGCGGCTGTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001142576

Insert Size: 1701 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142576.1](#)

RefSeq Size: 2526 bp

RefSeq ORF: 1701 bp

Locus ID: 3614

UniProt ID: [P20839](#)

Cytogenetics: 7q32.1

Protein Families:	Druggable Genome
Protein Pathways:	Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism
MW:	60.9 kDa
Gene Summary:	<p>The protein encoded by this gene acts as a homotetramer to regulate cell growth. The encoded protein is an enzyme that catalyzes the synthesis of xanthine monophosphate (XMP) from inosine-5'-monophosphate (IMP). This is the rate-limiting step in the de novo synthesis of guanine nucleotides. Defects in this gene are a cause of retinitis pigmentosa type 10 (RP10). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008]</p> <p>Transcript Variant: This variant (4) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (dz) has the same N- and C-termini but is shorter compared to isoform a.</p>