

Product datasheet for **SC326075**

TMEM151B (NM_001137560) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM151B (NM_001137560) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM151B
Synonyms:	bA444E17.5; C6orf137; TMEM193
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326075 representing NM_001137560.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCCCCTGGCTCGGCCGCGGGAGAGAGCGCCGCCGCGCGCGCGCGGTGGCGCCCCGGG
GTCTCGGAGGAGCTACGGCGCGCGCGCAGCGCGCGCGGACGAGGGCCCCCGCCGAGAGGAGCAG
CGTCCCATCCAGCCCTTTTACCAAGTCCCTCTGCCGTGAGTCCCACTGGAAAGTGCCTCCTGCTCG
CTGCTCATGTACGGCTGCCTGGGGCAGTGGCCTGGTGCCACGTACACAGTGACGCGCCTCACCTTC
AGCAGCGCTACCGGGCAACAGCCTCATGTACCATGACAGCCCCCTGCTCCAACGGCTATGTCTACATC
CCCCTGGCCTTCTGCTCATGTTGTACGCCGTCTACCTGGTGGAGTGTGGCACTGCCAAGCCGCCAT
GAGCTGCAGCACCCTGTTGATGTGAGCAGTGTGCGGGAACGTGTGGGCCCATGCAGCAAGCCACGCC
TGCATCTGGTGAAGGCCATCAGTACCACTATGTCCGCCGACCCGCCAGGTACACAGATACCGCAAT
GGAGACGCTATACCAACCCAGGTCTACCACGAACCGTCAACACGCACGTGGCGGAGGTGAGTTC
GACTACGCGCGCTGCGGCGTTCCGACGTGTCAAGACGTGGTGGGGCTGGAGGGCGCGCGGCCACG
CGGCTGCGCTTACCAAGTGCTTCAGTTTCGCCAGCGTGGAGGCCGAGAACGCGTACCTGTGCCAGCGC
GCGCGCTTCTTCGAGAGAACGAGGGCCTAGACGACTACATGGAGGCACGCGAGGGCATGCACCTCAAG
AACGTGGACTTCCGTGAGTTCATGGTGGCCTTCCCGGACCCGGCCCGGCCGCCCTGGTACGCTGCTCG
TCGGCCTTCTGGGCCGCGCGCTGCTCACGCTGTCTGGCCGCTGCGAGTGTGGCCGAGTACCGCACG
GCCTACGCGCACTACCAGTGGAGAAGCTATTTGGCCTGGAGGGCCCGGGCTCGGCCAGCAGCGCAGGC
GGTGGCCTCAGCCCCAGCGATGAGCTGTGCCCCGCTCACCACCGCTGCCGCGGGTCAACACAGTA
GACAGCACGGAGCTCAGTGGACATCCGCTCCAACAGCAGCTGGTGCCAGCTACTCTGAGGCGGTG
CTCATGGACCTGGCGGGCTCGGGACGCGCTGCGCGGGGAGGCGCGGCTACGCGCCCTCGTGCCGC
TACGGTGGGTAGCGGGCCCGCGCGCGCGGCGTGGCTCCCTACCGGCGCAGCTGCGAGCACTGCCAG
CGCGCGCTCAGCAGCTCGTCTATCTTCTCGCGCAGCGCCCTAAGCATCTGCGCCAGCCCGCGGGCGGC
CCGGGGCCCGGTGGGGCGCGGGCTGCGGGGAGCGCTTCTGCTGGCCGTCTCTACGGCTCCCGG
CGCAGCTGCCTGTGGCGAGCCGACGCGGGAGCGTCAACGAGGCCAGTGGCCACGAGCAGACGCGG
CTGTCCAGCCAGGCCAGCATGGGGGACGACGAGGACGACGAGGAGGAGGCGGGCCGCCGCCGCC
TACCACGACGCCCTCTACTTTCCGGTCTCATCGTCCACCGCAGGAGGGGTGTCTGGGCCACAGCCAC
CGGCCGCTGCACCGCCACGGCTCTGCGTAGAGACCTACTGTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTAA
  
```

Restriction Sites: Sgfl-RsrII

ACCN: NM_001137560

Insert Size: 1701 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001137560.1

RefSeq Size: 4637 bp

RefSeq ORF: 1701 bp

Locus ID: 441151

UniProt ID: Q8IW70

Cytogenetics: 6p21.1

Protein Families: Transmembrane

MW: 61.5 kDa