

## Product datasheet for **SC326079**

### **ABLIM2 (NM\_001130085) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ABLIM2 (NM_001130085) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABLIM2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC326079 representing NM\_001130085.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTGCAGTGTCCGAGCCCCAGGCTGCTCCAGCCCCGCTGGAGAAGTCGCCAGCACGGCGATCCTG
TGCAACACCGTGTGGGAATGTGTGCAAGGGCGAGGTGCTGCGGGTGCAGGACAAGTACTTCCACATCAAG
TGCTTCGTCTGTAAGCATGTGGCTGCGACCTGGCCGAGGGCGGCTTCTTCGTGCGGCAGGGCGAGTAC
ATCTGCACGCTGGACTACCAGAGGCTCTACGGACCCCGCTGCTTCAGCTGCGACCAGTTCATTGAGGGT
GAGGTGGTGTGCGGCTGGGCAAGACCTACCACCCCGACTGCTTCGTGTGTGCCGTCTGCCGGCTGCC
TCCCCCCCCGGGACCGAGTGACCTTCAACGGGAAGGAATGCATGTGCCAGAAGTGTTCCTGCCGTA
TCGGTGGGCGAGCAGCGCACCTGTCCCAGGGCTCCGAAGTGTGGGGCTGCGGCACAGAAATCAAG
AATGGCCAGGCCCTGGTAGCCTTGACAAGCACTGGCACTTGGGTGTTTTAAGTGAAGAGCTGTGGG
AAGCTCCTGAATGCCGAGTACATCAGCAAGGATGGGTGCCCTACTGCGAAGCTGACTATCACGCCAAG
TTCGGCATCCGCTGTGACAGCTGTGAGAAATACATCACGGGGCGGTGCTGGAGGCCGGAGAGAAGCAC
TACCACCTTCTCGCGCTATGTGTCAAGTGCAGCCAGATGTTTGCAGAAGGCCAAGAGATGTATCTT
CAAGGTTCTCCATCTGGCATCCGGCGTGTGACAAGCAGCCAGAAGTGAAGACAGAAACAGGAAACC
AGAATTCCTCAGAGAGCATCATTTCTGTCCCTGCTTCCAGCACCTCAGGGTCTCCGAGCCGTGTGATT
TATGCCAAGCTTGGTGGTGAATCCTGGACTACAGGGACTTGGCAGCCCTTCTAAAAGTAAAGCCATC
TATGACATCGACCGCCCCGACATGATCTCCTACTACCCTACATCAGCCACTCTGCAGGGGACAGGCAG
AGCTACGGCGAGGGGATCAGGATGACCGGTCTACAAGCAGTGTGCGACCTCCAGCCAAAGCTCCACT
GGTTCGGTTAGCCTCGGGCGTACACTCCGACCTCACGGTACCACAGCACTACAGCCGTCACGGCAGT
GAAAGTGGCCGGAGCACCCCGACCTCTCCGTCTCTGACAGCAAGCCGCCCCCTCCACCTACCAG
CAGGCACCTCGCCACTTCCACGTCCAGACTGGCGTAAAAGATAACATCTATAGGAAACCCCTATC
TACAGACAGCATGCTGCCAGCGATCGGATGGGGAGGATGGAAGCTTGGACCAGGATAACAGGAAGAAG
AGCAGCTGGTGTGCTCAAGGGGGATGCAGACACAAGGACCAATTCTCCAGACCTGGACACCCAGTCC
TTGTCCCACAGCAGCGGGACCGACAGAGACCTCTCCAAAGGATGGCAGGGGACAGCTTTCACTCACA
TACAAGATCTATCCGTATGACTCCCTCATCGTCACAAACCGAATTCGCGTGAAGTCCCAAAGACGTG
GACCGGACGAGACTGGAGAGACTTGTGCGCCGAGGAGTCCAGGAAGTGTGGGATGAGCATCGAG
GAGTTTACCAGCCTGGCCCTCTGGAAGAGGAATGACCTTAAGAAGAAAGCCCTTTTGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001130085
<b>Insert Size:</b>	1719 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001130085.1</a></u>
<b>RefSeq Size:</b>	3555 bp
<b>RefSeq ORF:</b>	1719 bp
<b>Locus ID:</b>	84448
<b>UniProt ID:</b>	<u><a href="#">Q6H8Q1</a></u>
<b>Cytogenetics:</b>	4p16.1
<b>Protein Pathways:</b>	Axon guidance
<b>MW:</b>	63.6 kDa
<b>Gene Summary:</b>	<p>May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon and uses an alternate in-frame splice site in the mid-coding region, and lacks a segment in the 3' coding region, compared to variant 1. The resulting isoform (3) is shorter than isoform 1.</p>