

Product datasheet for **SC326080**

Epsin 1 (EPN1) (NM_001130072) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Epsin 1 (EPN1) (NM_001130072) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPN1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326080 representing NM_001130072.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCGACCTCGTCCTTGAGGCGCCAGATGAAGAACATCGTCCACAACACTACTCAGAGCGGAGATCAAG
GTTTCGAGAGGCCACGAGCAATGACCCCTGGGGCCCATCCAGCTCCCTCATGTCTCAGAGATTGCCGACCTC
ACCTACAACGTTGTCGCCTTCTCGGAGATCATGAGCATGATCTGGAAGCGGCTCAATGACCATGGCAAG
AACTGGCGTCACGTTTACAAGGCCATGACGCTGATGGAGTACCTCATCAAGACCGGCTCGGAGCGCGTG
TCGCAGCAGTGCAAGGAGAACATGTACGCCGTGCAGACGCTGAAGGACTTCCAGTACGTGGACCGCGAC
GGCAAGGACCAGGGCGTGAACGTGCGTGAGAAAGCTAAGCAGCTGGTGGCCCTGCTGCGCGACGAGGAC
CGGCTGCGGGAAGAGCGGGCGCACGCGCTCAAGACCAAGGAAAAGCTGGCACAGACCGCCACGGCCTCA
TCAGCAGCTGTGGGCTCAGGCCCCCTCCCGAGGCGGAGCAGGCGTGGCCGAGAGCAGCGGGGAGGAG
GAGCTGCAGCTCCAGCTGGCCCTGGCCATGAGCAAGGAGGAGGCCAGCAGCCCGTCTGCGGCCCC
GAGGACGACGCCAGCTCCAGCTGGCCCTTAGTTTGGCCGAGAAGAGCATGATAAGGAGGAGCGGATC
CGTCGCGGGGATGACCTGCGGCTGCAGATGGCAATCGAGGAGAGCAAGAGGGAGACTGGGGCAAGGAG
GAGTCGTCCTCATGGACCTTGCTGACGCTTTCACGGCCCGACTCTGCCCCGACCACAGACCCCTGG
GGGGGCCAGACCCATGGCTGCTGCCGTCCCCACGGCTGCCCCACCTCGGACCCCTGGGGCGGCCCC
CCTGTCCCTCCAGCTGCTGATCCCTGGGAGGTCCAGCCCCACGCCGGCCTCTGGGGACCCCTGGAGG
CCTGCTGCCCCTGCAAGACCCCTCAGTTGACCCTTGGGGTGGGACCCAGCCCTGCAGCTGGGGAGGGG
CCCACGCTGATCCATGGGAAGTTCCGATGGTGGGTCCCGTCAGTGGGCCCTCAGCCTCCGATCCC
TGGACACCGGCCCGGCCCTTCTCAGATCCCTGGGAGGGTCACTGCCAAGCCCAGCACCATGGCACA
ACAGCAGCCGGGGATTTCGACACGGAGCCCGACGAGTTCTTGACTTTGACCGACTCCGACGCGCACTG
CCGACCTCCGGGAGCAGCGCAGGAGAGCTGGAGCTGCTGGCAGGAGAGGTGCCGGCCGAAGCCCTGGG
GCGTTTGACATGAGTGGGTGAGGGATCTCTGGCTGAGGCTGTGGGAGCCCCACCTGCAGCCACA
CCAACCTCCACGCCCCCACCGGAAGACGCGGAGTCAATTCCTGGGGCCCAATGCAGCCCTCGTCGAC
CTGGACTCGTGGTGAGCCGGCCGGGCCACGCCGCTGGAGCCAAGGCCTCCAACCCCTTCTGCCA
GGCGGAGGCCAGCCACTGGCCCTTCCGTACCAACCCCTTCCAGCCCGCCTCCCGGACGCTCACC
CTGAACCAGCTCCGTCTCAGTCTGTGCCTCCCGTCCCTGGAGCGCCACCCACGTACATCTCTCCCTT
GGCGGGGCCCTGGCTGCCCCATGATGCCCCGGGCCCCCGGCCCAACACTAATCCCTTCTCT
CTATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001130072

Insert Size: 1731 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130072.1](#)

RefSeq Size: 2501 bp

RefSeq ORF: 1731 bp

Locus ID: 29924

UniProt ID: [Q9Y6I3](#)

Cytogenetics: 19q13.42

Protein Pathways: Endocytosis

MW: 60.3 kDa

Gene Summary: This gene encodes a member of the epsin protein family. The encoded protein binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Loss of function of this gene is associated with reduced tumor growth and progression in certain cancer types. [provided by RefSeq, Mar 2016]

Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and includes an additional internal exon, compared to variant 1. The resulting isoform (b) is shorter at the N-terminus, compared to isoform a.