

Product datasheet for **SC326084**

ACOXL (NM_001142807) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACOXL (NM_001142807) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACOXL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326084 representing NM_001142807.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAGAGCTTTGACAGTTCAGAGAGTGAAGTTTGCCATGGACCTGCCTCTGTAAAAACGTGCAGGTCAG
GATCTGGCAGAGAAAACAAAGAATTTTGTACGCCGAAGCCTTGTCATAGGAGAAGTCCTCTCCATGGCG
GACATGGCCACAGGAGTGAAGTGCGAATAATTTATTGGCTATTTGGTGGTGCTATCAGGAATCTCGGA
AGCCCTGAACATGTTACTAAGTGGTTTCAGCCACTCCAGGAGCAGAAATACACTGGGATGTTTGCAATG
ACCGAGAGGGGCCATGGGAGCAACGCGAGAGGGATCCAGACCGAAGCCACCTTTGACCTCTCTGCCAG
GAGTTTGTAAATTGACACGCCGTGTGAAATGCGGAGAAGATGTATATTGAAATGCCATGTACGGGAAT
TATGCAGCTGTCTTTGCCAGCTCATCATAGATGGAAGATCTCAAGGGCCCCACTGTTTCATCGTTCCT
GTCCGGGATGAAACGGAAGCTTGTACCCAGGAGTCACAGCTATTGATATGATGTACAAGGAGGGTCTG
CATGGTGTGGACAATGGGATATTAATTTTACAAGGTTTCGGATACCCAGGAGAACCTGCTGGATAAG
TTTGGTCCGTGGCTCCAGATGGACAGTACCATTTCGCTATTAGGAACAAGAGTGCAAGATTCAATGCC
ATGCTGGCAGCACTGACCCCTTCGAGATTAGCTGTGGCTTCCAAGCTATGGGTGCCATGAAGCTTGGG
TTGACGATAGCCATTTCGTATAGCCACAGCCGGAGGCAGTTTGGGCCAAAACCAAGGAAGAGGTGAAG
ATCATTGAGCACCACAAACAGACCCCTGCCGCTGATGCCCCACCTGGCCACAGCCTTGGCCCTGACCTTC
GTCAGCAGGTATGCTGGGGCCCTCCTGGATGAGGATGTCTTCAGGGAAGGAGCTGGTCAACAGTCGC
TCGCTGCAGGCTCTGGTGGCGGGCTGAAGGCCTACAGCACCTGGGAGAACATCCGCTGCCTGCAGGAC
TGCCGCGAGTGCCTGGAGGCATGGTTGTGGGGCGGGAAGTCTGGCCCAATACACCAACAGTATGAA
GAAAAACCACTCTTTGGCCTGCTCCAAACTGGGCTGAATCTGTGGGGGACAAGCTGAGAACCAGTTTC
CTGGCATTAAACATGGACACAGTTGATGATCTCGCCTTTCTGTTGAAAGCAGTGAAATTTCTGTAAGG
GTTCTTCAGCGGGGTTTGGTGGCCAGAATTTATTATAAGGTAAAGACCAAGAAAGGAGGATTTTTCCAT
GCCTGGAAGTCTGTCTGCACCAGTGGCTTCTGTCCCTGGCACACACTCACCAGTTACCTTAGAG
CAGTTCTCCCTAGCAGTGAAGAGCTGTCTGACCAAGAGGACCAGACTTTGTTAATGAAGTTTGTCTG
TTGTATGGAACCAAGCTGGTGTTCAGGAGCGGCCTGGTATTTAGAACATAAATACTTGACTCCCATG
GCCAGCACGAGGATCAGGAATCAGTTGCTGGATTTGTGCGACTCGGTGAAGGATGATGCCCGAGGGTG
ATCTCGACCTTTAACATTCCACACACCTACCTCCACGCACCAATCGCCGAATCTCCAACCCGCGGGCC
GCGTGGGCTTTCTACCCTGCACCGCTGCAGCCGCGGCCACGGGAAGAGGCGCGCTCCCGCGGCCCAAG
CTGGGAGCCAAGCTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001142807

Insert Size: 1743 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142807.1](#)

RefSeq Size: 2373 bp

RefSeq ORF: 1743 bp

Locus ID: 55289

UniProt ID: [Q9NUZ1](#)

Cytogenetics: 2q13

MW: 65.5 kDa