

Product datasheet for **SC326085**

TTC39A (NM_001144832) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39A (NM_001144832) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39A
Synonyms:	C1orf34; DEME-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326085 representing NM_001144832.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCTTTCTGACCACCTGGAAGGAGGCGCCCTCCGGGAAAAGCGACCGCAGGACTCCTGAGAGCAGC
CTCCATGAGGCCCTGGACCAGTGCATGACCGCCCTGGACCTTCTCCTACCAACCAGTTCTCAGAAGCA
CTCAGCTACCTCAAGCCCAGAACCAAGGAAAGCATGTACCACTCACTGACATATGCCACCATCCTGGAG
ATGCAGGCCATGATGACCTTTGACCCTCAGGACATCCTGCTTGCCGGCAACATGATGAAGGAGGCACAG
ATGCTGTGTCAGAGGCACCGGAGGAAGTCTTCTGTAACAGATTCTTCAGCAGCCTGGTGAACCGCCCC
ACGCTGGGCCAATTCACTGAAGAGGAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGA
GCAGCCCTGACCTTCTGCAGGACGAGAACATGGTGAAGTTCATCAAAGGCGGCATCAAAGTTCGAAAC
AGCTACCAGACCTACAAGGAGCTGGACAGCCTTGTTCAGTCTCACAACTGCAAGGGTGAGAACCAC
CCGCACCTTGAAGGAGGAGTGAAGCTTGGTGTAGGGGCCCTCAACCTGACTGTCCATGCTTCTACT
AGGATCCTGAGGCTGTTGGAGTTTGTGGGTTTTTCAGGAAACAAGGACTATGGGCTGCTGCAGCTGGAG
GAGGGAGCGTCAGGGCACAGCTTCCGCTCTGTGCTCTGTGTCATGCTCCTGCTGTGCTACACACCTTC
CTCACCTTCGTGCTCGGTACTGGGAACGTCAACATCGAGGAGGCCGAGAAGCTCTTGAAGCCCTACCTG
AACCGGTACCCCTAAGGGTGCCATCTTCTGTTCTTTCAGGGGAGGATTGAAGTCATTAAGGCAACATT
GATGCAGCCATCCGGCGTTTCGAGGAGTGTGTGAGGCCAGCAGCACTGGAAGCAGTTCACCACATG
TGCTACTGGGAGCTGATGTGGTGCTTACCTACAAGGGCCAGTGGAAAGATGTCCTACTTCTACGCCGAC
CTGCTCAGCAAGGAGAACTGCTGGTCCAAGGCCACCTACATTTACATGAAGGCCGCTACCTCAGCATG
TTTGGGAAGGAGGACCACAAGCCGTTCCGGGACGACGAAGTGAATTATTTGAGCTGTGCCAGGCCTG
AAGCTCAAGATTGCTGGGAAATCTCTACCACAGAGAAGTTTGCCATCCGGAAGTCCCGGCCACTTTC
TCTCCAACCCTATCTCGCTGCCAGTGCCTGCTCTGGAATGATGTACATCTGGAACGGCTACGCCGTG
ATTGGGAAGCAGCCGAAACTCACGGATGGGATACTTGAGATTACTACTAAGGCTGAAGAGATGCTGGAG
AAAGGCCCAGAGAACGAGTACTCAGTGGATGACGAGTGTGGTGAATTTGTTGAAAGGCCCTGTGCTG
AAATACCTGGGCCGTGTCCAGGAGGCCGAGGAGAAATTTAGGAGCATCTCTGCCAATGAAAAGAAGATT
AAATATGACCACTACTTGATCCCAAACGCCCTGCTGGAGCTGGCCCTGCTGCTTATGGAGCAAGACAGA
AACGAAGAGGCCATCAAACCTTTGGAATCTGCCAAGCAAACACTACAAGAATTACTCCATGGAGTCAAGG
ACACACTTTCGAATCCAGGCAGCCACACTCCAAGCCAAGTCTTCCCTAGAGAACAGCAGCAGATCCATG
GTCTCATCAGTGTCTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001144832

Insert Size: 1746 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001144832.2](#)

RefSeq Size: 2736 bp

RefSeq ORF: 1746 bp

Locus ID: 22996

UniProt ID: [Q5SRH9](#)

Cytogenetics: 1p32.3

MW: 66.6 kDa