

Product datasheet for **SC326087**

TTC39C (NM_001135993) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39C (NM_001135993) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39C
Synonyms:	C18orf17; HsT2697
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326087 representing NM_001135993.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCGGCTCGGAGCAGCAGCGGCCGCGGCGGGACGACGGAGACTCGGACGCGGCAGCGCGGGC
GCGGCGCCCTGCAGGACGCGGAGCTGGCCCTGGCCGGCATCAACATGCTGCTCAACAACGGCTTCAGG
GAGTCGGACCGAGCTTTTCAAACAATACAGAAATCATAGCCCACTAATGAGTTTTGGAGCCAGCTTTGTC
AGTTTTTTGAATGCCATGATGACATTTGAGGAAGAAAAATGCAGTTGGCATGTGATGACTTAAAAACC
ACAGAAAACTGTGTAAAGTGAAGAGGCTGGAGTAATTGAAACAATCAAGAATAAAATTAAGAAGAAC
GTTGATGTCCGAAAAATCCGCCCCCTCTATGGTTGATCGGCTTCAGAGGCAGATAATCATAGCTGACTGC
CAGGTTTACCTGGCTGTGCTTTTCAATTTGTAACAAGAATTGTCAGCTTATCAAGGTGGGTGGATC
CTTAGGAAAGCCTGGAAGATTTACAATAAATGCTATCTGGACATCAATGCCCTTCAGGAGCTGTATCAG
AAGAAGCTAACTGAAGAGTCCTTGACTTCTGATGCTGCAATGATAATCACATTGTGGCTGAAGGGGTG
TCTGAGGAGTCTCTGAACAGACTGAAAGGTGCTGTTAGCTTTGGATATGGCCTTTTTCACCTTTGCATA
TCCATGGTGCCCCCAAACTGCTCAAAATCATCAACCTGCTGGGTTTTCTGGAGACCGCTACAGGGG
CTTTCTTCACTGATGATGCAAGCGAAAGTAAGGACATGAAGGCCCTTTAGCTACATTAGCTCTGCTC
TGGTATCATACTGTAGTCCGCCGTTTTTGCCTTGGATGGCAGTGATAACAAGGCAGGCTGGATGAA
GCTAAGGAAATTCCTTAAAAAAGAAGCTGCTTATCCAAATCTTCCCTCTTTATGTTTTCAAGGGA
CGGATACAACGACTAGAGTGTCAATCAACAGTGCCTTGACATCTTCCACACTGCTTTGGAAGTTGCA
GTAGACCAGAGAGAAATCAACATGTCTGTCTGTATGAAATTGGTTGGTGCAGCATGATAGAGCTCAAT
TTCAAGGATGCATTTGATTCCTTTGAGAGGCTAAAAATGAGTCCAGGTGGTCCAGTGCTATTATGCC
TACTTGACTGCAGTTTGTGAGGAGCCACTGGTGTGATGTGGATGGGGCACAGATTGCTTTAAAGAAGTT
CAGAAACTCTTCAAAAGGAAAAACAATCAGATTGAACAGTTCTCGGTGAAAAAGGCAGAGCGATTTCGG
AAGCAAAACCCCAACCAAGCGCTCTGTGTGTTGGCGTCTATTGAAGTGTGACTTGTGGAAGCTCTT
CCAAACTGTTCTTCCCAACCTGCAGAGGATGAGTCAAGCTTGCCATGAAGTGGATGACTCATCTGTT
GTTGGATTAAAGTATTTGCTTCTTGGTGCCATACACAAATGTCTAGGAACTCAGAAGATGCTGTTGAG
TACTTCCAGCGAGCTGTTAAAGATGAATTGTGTCGTGAGAATAATTTATATGTTGAGCCGTATGCCTGT
TATGAAGTGGTGTCTTCTATTAGACAAACAGAGACTGTAGGAAGAGGCAGAGCTCTACTTCTCAA
GCAAGGAGGATTCTCTGGCTACGACTTTGAAAACAGATTGCATGTCCGCATCCATGCTGCTCTGGCC
TCTCTGAGGGAATTGGTTCCTCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001135993

Insert Size: 1752 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135993.1](#)

RefSeq Size: 5253 bp

RefSeq ORF: 1752 bp

Locus ID: 125488

UniProt ID: [Q8N584](#)

Cytogenetics: 18q11.2

MW: 65.9 kDa