

Product datasheet for **SC326106**

ZFY (NM_001145276) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZFY (NM_001145276) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZFY
Synonyms:	ZNF911
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326106 representing NM_001145276.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATGAAGATGAATTTGAATTGCAGCCACAAGGCCAACTCATTTTTTATGGAATAGTGGATGAT
GCTGGCAAAATAGAACATGATGGTTCCACTGGAGTGACCATCGATGCAGAATCAGAAATGGATCCTTGT
AAAGTGGATAGCACTTGTCTGAAGTCATCAAGGTGTACATTTTAAAGCTGACCCCTGGAGAAGATGAC
TTAGGTGGAAGTGTAGACATTGTGGAGAGTGAACCTGAAAATGATCATGGAGTTGAACTACTTGATCAG
AACAGCAGTATTCGTGTTCCAGGGAAAAGATGGTTTATGACTGTCAATGACTCTCAACAAGAAGAT
GAAGATTTAAATGTTGCTGAAATTGCTGATGAAGTTTATGGAAGTGATCGTAGGAGAGGAGGATGCT
GCTGTTGCAGCAGCAGCAGCTGTGTCATGAGCAGCAAAATTGATGAGGATGAAATGAAAACCTTCGTA
CCAATTGCATGGGCAGCAGCTTATGGTAATAATTCTGATGGAATTGAAAACCGGAATGGCACTGCAAGT
GCCCTCTTGCACATAGATGAGTCTGCTGGCCTTGGCAGACTGGCTAAACAGAAACCAAGAAAAAGAGA
AGACCTGATTCCAGGCAGTACCAACAGCAATAATTATTGGCCCTGATGGTCATCCTTTGACTGTCTAT
CCTTGCATGATTTGTGGGAAGAAGTTTAAAGTCGAGGGGTTTTTAAAAAGACACATGAAAAACCATCCT
GAACACCTTGCCAAGAAGAAGTACCACTGTACTGACTGTGATTACACTACCAATAAGAAGATAAGTTTA
CATAACCACCTGGAGAGCCACAAGCTGACCAGCAAGGCAGAGAAGGCCATTGAATGTGATGAGTGTGGG
AAGCATTTTTCTCATGCAGGGGCTTTGTTTACTCACAATGGTGCATAAGGAAAAAGGGGCCAACAAA
ATGCACAAGTGTAATTTCTGTGAATATGAGACAGCTGAACAGGGGTTATTGAATCGCCACCTCTTGGCA
GTCCACAGCAAGAAGTTTCTCATATTTGTGTGGAGTGTGGTAAAGGTTTCCGACACCCGTCGGAAGT
AGAAAGCAGATGCGAATCCATACCGGCGAGAAGCCATACCAATGCCAGTACTGTGAATATAGGTCTGCA
GACTCTTCTAAGTTGAAAACACATATAAAAAACAAAGCATAGTAAAGAGATGCCATTCAAGTGTGACATT
TGTCTTCTGACTTTCTCAGATACCAAGAAGTGCAGCAACATACTCTTGTCCACCAAGAAGCAAAACA
CATCAGTGTGTTGCATTGCGACCAAGAGTTCAAAGTCAAGTATTGAAACGACATGTAATTTAGTT
CATACGAAAGACTATCCTCATAAGTGTGAGATGTGCGAGAAAGGCTTTCACAGGCCTTCAGAACTTAAG
AAACATGTGGCTGTCCACAAAGGTAAAAAATGCACCAATGTAGACATTGTGACTTTAAGATTGCAGAC
CCATTTGTTCTAAGTCGCCATATTCTCTCAGTTACACAAAGGATCTTCCATTTAGGTGTAAGAGATGT
AGAAAGGGATTTAGGCAACAAAATGAGCTTAAAAAGCATATGAAGACACACAGTGGCAGGAAAGTATAT
CAGTGTGAGTACTGTGAGTATAGCACTACAGATGCCTCAGGCTTTAAACGGCACGTTATTTCCATTCAT
ACAAAAGACTATCCTCATCGGTGTGAGTACTGCAAGAAAGGCTTCCGAAGACCTTCAGAAAAGAACCA
CACATAATGAGACACCATAAAGAAGTTGGTCTGCCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001145276

Insert Size: 1833 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145276.1</u>
RefSeq Size:	4769 bp
RefSeq ORF:	1833 bp
Locus ID:	7544
UniProt ID:	<u>P08048</u>
Cytogenetics:	Yp11.2
Protein Families:	Transcription Factors
MW:	69.9 kDa
Gene Summary:	This gene encodes a zinc finger-containing protein that may function as a transcription factor. This gene was once a candidate gene for the testis-determining factor (TDF) and was erroneously referred to as TDF. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (3) includes an alternate 5'UTR exon and omits an in-frame coding exon, compared to variant 1, resulting in a shorter protein (isoform 3).