

Product datasheet for **SC326110**

Ribophorin II (RPN2) (NM_001135771) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ribophorin II (RPN2) (NM_001135771) Human Untagged Clone
Tag:	Tag Free
Symbol:	RPN2
Synonyms:	RIBIIR; RPN-II; RPNII; SWP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326110 representing NM_001135771.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGCCGCGGGTTCAAGCACTGTCTTCTGTTGGCCCTGACAATCATAGCCAGCACCTGGGCTCTG
ACGCCCACTCACTACCTACCAAGCATGACGTGGAGAGACTAAAAGCCTCGCTGGATCGCCCTTTCACA
AATTTGGAATCTGCCTTCTACTCCATCGTGGGACTCAGCAGCCTTGGTGCTCAGGTGCCAGATGCAAAG
ATCTCTATTTCAAATGAGACCAAAGATCTGCTTCTGGCAGCTGTCAAGTGGACTCATCTGTTACCCAG
ATCTACCATGCAGTTGCAGCTCTAAGTGGCTTTGGCCTTCCCTTGGCATCCCAAGAAGCACTCAGTGCC
CTTACTGCTCGTCTCAGCAAGGAGGAGACTGTGCTGGCAACAGTCCAGGCTCTGCAGACAGCATCCAC
CTGTCCCAGCAGGCTGACCTGAGGAGCATCGTGGAGGAGATTGAGGACCTTGTGCTCGCCTGGATGAA
CTCGGGGGCGTGTATCTCCAGTTTGAAGAAGGACTGGAAACAACAGCGTTATTTGTGGCTGCCACCTAC
AAGCTCATGGATCATGTGGGACTGAGCCATCCATTAAGGAGGATCAGGTATCCAGCTGATGAACGCG
ATCTTCAGCAAGAAGAATTTGAGTCCCTCTCCGAAGCCTTCAGCGTGGCCTCTGCAGCTGCTGTGCTC
TCGCATAATCGCTACCAGTGCCAGTTGTGGTTGTGCCTGAGGGCTCTGCTTCCGACACTCATGAACAG
GCTATCTTGCGGTTGCAAGTCAACAATGTTCTGTCTCAGCCTCTGACTCAGGCCACTGTTAAACTAGAA
CATGCTAAATCTGTTGCTTCCAGAGCCACTGTCCTCCAGAAGACATCCTTCAACCCTGTAGGGGATGTT
TTTGAACATAAATTTATGAACGTCAAATTTTCCAGTGGTTATTATGACTTCCCTTGTGCAAGTTGAAGGT
GACAACCGGTATATTGCAAATACCGTAGAGCTCAGAGTCAAGATCTCCACTGAAGTTGGCATCACAAT
GTTGATCTTCCACCGTGGATAAGGATCAGAGCATTGCACCCAAAACACCCGGGTGACATACCCAGCC
AAAGCCAAGGGCACATTCATCGCAGACAGCCACCAGAATTCGCCTTGTCTTCCAGCTGGTAGATGTG
AACACTGGTGTGAACACTCCTCACCAGACATTTGTCCGACTCCATAACCAGAAGACTGGCCAGGAA
GTGGTGTGTTGTTGCCGAGCCAGACAACAAGAACGTGTACAAGTTTGAACCTGGATACCTCTGAAAGAAG
ATTGAATTTGACTCTGCCTCTGGCACCTACACTCTCTACTTAATCATTGGAGATGCCACTTTGAAGAAC
CCAATCTCTGGAATGTGGCTGATGTGGTCAAGTTCCTGAGGAAGAAGCTCCCTCGACTGTCTTG
TCCCAGAACCTTTTCACTCCAAAACAGGAAATTCAGCACCTGTTCCGCGAGCCTGAGAAGAGGCCCCCC
ACCGTGGTGTCCAATACATTCCTGCCCCTGATCCTCTCGCCGTTGCTTCTGCTCTCGCTCTGTGGATC
CGGATTGGTGCCAATGTCTCAACTTCACTTTTGTCTCCTAGCACGATTATATTTACCTGGGACATGCT
GCTATGCTGGGACTCATGTATGTCTACTGGACTCAGCTCAACATGTTCCAGACCTTGAAGTACCTGGCC
ATCCTGGGCACTGTGACGTTTCTGGCTGGCAATCGGATGCTGGCCAGCAGGCAGTCAAGAGGATCGT
GCAGAGCAGAGCAGTAGATTGGCAAAATATAGGACACTACGAACAGCACATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001135771

Insert Size: 1848 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135771.2](#)

RefSeq Size: 2497 bp

RefSeq ORF: 1848 bp

Locus ID: 6185

UniProt ID: [P04844](#)

Cytogenetics: 20q11.23

Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis
MW:	67.7 kDa
Gene Summary:	<p>This gene encodes a type I integral membrane protein found only in the rough endoplasmic reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus motif of nascent polypeptide chains. This protein is similar in sequence to the yeast oligosaccharyl transferase subunit SWP1. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (2) encodes isoform 2. Variants 2 and 3 are the same length and encode isoforms (2 and 3, respectively) of the same length but their variant and encoded isoform sequences are different.</p>