

Product datasheet for **SC326128**

NGLY1 (NM_001145293) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NGLY1 (NM_001145293) Human Untagged Clone
Tag:	Tag Free
Symbol:	NGLY1
Synonyms:	CDDG; CDG1V; PNG-1; PNG1; PNGase
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326128 representing NM_001145293.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGCGCGGCATTGGGCAGCTCCTCAGGCTCGGCGTCCCCGGCCGTGGCTGAGCTCTGCCAGAAC
ACCCCGGAGACCTTTTTGGAGGCTCCAAGCTGCTGCTACCTATGCTGACAACATCCTCAGAAACCTT
AATGATGAAAAATATAGATCCATCCGGATTGGAAACACAGCCTTTTCTACTAGACTCTTGCTGTGAGA
GGAGCTGTTGAATGTTTATTTGAAATGGCTTTGAAGAGGGAGAAACACATCTCATCTTTCTAAAAAA
GCTTCAGTGGAGCAGCTGCAAAAAATTCGTGACCTGATTGCCATAGAGAGAAGTAGCAGACTGGATGGC
TCAAAATAAGAGCCACAAAGTAAAGTCACTCAGCAACCTGCAGCCAGTACCCAGCTTCTACAACACCA
TCTTCAAATCCCAGTGGTTAAACCAGCACACAAGGAACCGTCAAGGGCAGTCATCAGATCCACCATCT
GCTTCAACGGTTGCTGCTGACTCAGCCATTCTAGAAGTCTTCAGTCCAACATTCAGCATGTGCTGGTC
TATGAAAATCCTGCTCTCAGGAGAAAGCGTTGGCTTGATTCCGGTCCAAGAACTAAAAAGGAAATCA
CAAGAAAAGTTATCGAGAGCTAGAAAATTGGATAAAGGTATCAATATAAGTGATGAGGATTTCTTTTG
CTGGAGCTTTTGACTGGTTAAGGAAGAATTTTTCACTGGGTGAATAACGTTTTGTGCAGCAAAATGT
GGTGGACAGACTAGGTCTAGAGATAGATCATTACTGCCAGTGATGATGAGCTGAAGTGGGTGCAAG
GAAGTGGAAAGATCACTACTGTGATGCCTGCCAGTTCAGCAATCGATTCCCAAGATATAATAACCCTGAG
AACTTTTGGAAACAAGATGTGGACGGTGTGGCAGTGGGCCAATTGTTTTACACTGTGCTGCCGAGCT
GTAGGGTTTGAAGCTCGCTATGTTTGGGATTACACAGAGTTACAGCGCACATTAAGCCTAAAACTGACC
ACGTTAAAGGAGATTGAAAAACATTTTTAAGGATTTCAAAGAGACAGAAATTAATACAGGTAGTTGAT
GTCACCTGGCGATATTCCTGCAACATGAAGAGGTGATTGCCAGAAGAACTAAGGTTAAAGAAGCATTAA
CTTCGAGACACTATTAATGGCTTAATAAGCAGAGGCAACTGTTTTTGTGAGAAAACAGAAGGAAAGAA
CTTCTCCAGAGGATAATTGTGGAGCTTGTGAATTTATATCTCCAAAACCCCTAAACCTGGAGAATTT
GGGGAAGAATATCTGGGTGAGTGGCTTGAGAGTAGCCCGAGGTGAAATGGGTCTACAGAGAAAAGAA
ACCTTGTTTATTCCTGTGAAAATGAGAAGATTTCTAAACAGCTCCACCTTTGTTACAATATTGTGAAA
GATCGTTATGTTGAGTTTCAAATAACAATCAAACCATTTCTGGATGGGAGAATGGCGTGTGAAAATG
GAATCTATATTCAGAAAAGTTGAAACAGACTGGCACATGGTATATTTGGCCGAAAGGAAGGATCATCT
TTTGCTTATATTTCTGGAAGTTTGAGTGTGGTCACTGGCCTAAAAGTAGATAGCATTCTATTAGA
ACAAGTAGTCAAACCTTTTCAGACTGGAACAGTAGAATGGAAATTCGATCTGATACAGCACAAAGTAGAA
CTGACAGGCGATAACAGTCTTCACTCCTATGCTGATTTTTCTGGTGCCACTGAAGTATTTTGGAAAGCA
GAATTAAGCAGAGGAGATGGTGTGTCGCTTGGCAACACACCCAGCTGTTTAGACAAAGCTTAAATGAC
CATGAAGAAAATTTGTTGGAGATAAATTAAAAATTCAGTGACCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001145293

Insert Size: 1911 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145293.1](#)

RefSeq Size: 2545 bp

RefSeq ORF: 1911 bp

Locus ID: 55768

UniProt ID: [Q96IV0](#)

Cytogenetics: 3p24.2

MW: 72.3 kDa

Gene Summary: This gene encodes an enzyme that catalyzes hydrolysis of an N(4)-(acetyl-beta-D-glucosaminy) asparagine residue to N-acetyl-beta-D-glucosaminyamine and a peptide containing an aspartate residue. The encoded enzyme may play a role in the proteasome-mediated degradation of misfolded glycoproteins. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2009]
Transcript Variant: This variant (2) includes an alternate exon that causes a frameshift and lacks an alternate exon which restores the reading frame, compared to variant 1. The resulting protein (isoform 2) is shorter compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.