

Product datasheet for **SC326142**

HBS1L (NM_001145158) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HBS1L (NM_001145158) Human Untagged Clone
Tag:	Tag Free
Symbol:	HBS1L
Synonyms:	EF-1a; eRF3c; ERFS; HBS1; HSPC276
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326142 representing NM_001145158.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCCCGGCATCGGAATGTTTCGAGGCTATAACTACGATGAAGATTTTGAAGATGATGATCTCTACGGC
CAGTCTGTAGAGGATGATTATTGTATTTTCGCCGTCAACAGCTCGTCTTTATTCATGCCTTGATCACATG
AGAGAGGTACTTGGAGATGCTGTGCCAGATGAAATATTAATTGAAGCAGTTCTGAAGAACAAGTTTGAT
GTGCAGAAGGCTTTGTCAGGGTCTGGAACAAGATAGAGTGCAGAGTTTGAAGGACAAGATGAGGCA
ACAGTATCTACAGGAAAGATAGCAAAAGGAAAACCAAGTAGATTCCCAGACATCGCGAAGTGAATCTGAA
ATTGTGCCAAAAGTTGCTAAAATGACTGTATCTGGAAAGAAGCAAATATGGGATTTGAAGTGCCTGGA
GTATCTTCTGAAGAAAATGGTCATAGTTTCCACACACCTCAAAAAGACCGCCATTGAAGATGCCATT
GCTTCTCCGATGTTCTTGAGACTGCTTCTAAATCTGCTAATCCACCCACACGATTCAAGCATCAGAA
GAGCAGAGTTCAACCCAGCACCGGTGAAAAGTCTGGCAAGCTGAGGCAGCAAAATAGATGTGAAGCGC
GAACTGGAGAAGCGGCAAGGAGGAAGCAGCTACTCAACTTAGTGGTATTGGTCAATGTTGATGCTGGG
AAAAGTACTCTGATGGGCCATATGCTTTATCTTCTGGGTAAATATAAACAAGAACTATGCATAAGTAT
GAACAGGAGTCTAAAAGGCTGGCAAAGCTTCGTTTGCATATGCATGGGTCTTGATGAAACTGGCGAA
GAAAGGGAAAGGGGAGTAACCATGGATGTTGGTATGACAAAGTTTGAAACCACAACCAAGTTATTACA
TTAATGGATGCTCCAGGCCATAAGGACTTCATTCCAAATATGATTACAGGAGCAGCCAGGCGGATGTA
GCTGTTTTAGTTGTAGATGCCAGCAGGGGAGAGTTTGAAGCTGGATTTGAGACTGGAGGACAAACACGA
GAGCATGGACTCTTGGTCCGTTCTCTGGGAGTGACGCAGCTTGCAGTTGCAGTTAATAAAAATGGATCAG
GTTAATTGGCAACAAGAAAGTTTTCAAGAGATTACTGAAAACCTGGGCACCTTCTTAAGCAAGCAGGT
TTCAGTCAAGTGAACCTCACAAAATGGTATAAAGGACTATGTTTATTAGAACAATGATTCTTTAAGCCT
CCCCAGCGATCTATTGACAAACCTTTTAGATTATGTGTGTCGATGTTTTCAAAGATCAAGGATCTGGA
TTTTGCATAACTGGTAAAATAGAAGCTGGTTATATCCAAACTGGTGACCGACTACTGGCAATGCCTCCT
AATGAAACTTGTACCGTGAAGGAATCACTCTGCATGATGAACCTGTCGACTGGGCGGCAGCAGGCGAT
CATGTTAGTCTTACTTTGGTTGGGATGGATATCATCAAAATCAATGTTGGCTGCATATTTTGTGGCCCC
AAAGTACCCATTAAGCTTGCACTCGTTTCAGAGCCCGAATCCTCATCTTAATATTGAAATTCCTATC
ACTAAAGGATTTCTGTGCTGTTACACTACCAAACCTGTCAGTGAACCCGCCGTTATTAACGATTGATT
AGTGTCTTAAACAAAAGCACGGGTGAAGTCACAAAGAAAAGCCTAAGTTTTTGACTAAAGGCCAGAA
GCATTGGTAGAGCTACAGACACAAAGACCAATAGCTCTTGAGCTATATAAAGACTTTAAAGAGCTGGGG
AGGTTTCATGCTACGTTACGGTGGTTCTACAATAGCTGCTGGTGGTGTCACTGAGATAAAAGAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001145158

Insert Size: 1929 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145158.1](#)

RefSeq Size: 7037 bp

RefSeq ORF: 1929 bp

Locus ID: 10767

UniProt ID: [Q9Y450](#)

Cytogenetics: 6q23.3

MW: 70.6 kDa

Gene Summary: This gene encodes a member of the GTP-binding elongation factor family. It is expressed in multiple tissues with the highest expression in heart and skeletal muscle. The intergenic region of this gene and the MYB gene has been identified to be a quantitative trait locus (QTL) controlling fetal hemoglobin level, and this region influences erythrocyte, platelet, and monocyte counts as well as erythrocyte volume and hemoglobin content. DNA polymorphisms at this region associate with fetal hemoglobin levels and pain crises in sickle cell disease. A single nucleotide polymorphism in exon 1 of this gene is significantly associated with severity in beta-thalassemia/Hemoglobin E. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, May 2009]

Transcript Variant: This variant (2) lacks an in-frame exon in the coding region, compared to variant 1. This results in a shorter protein (isoform 2) compared to isoform 1. **Sequence Note:** This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.