

Product datasheet for **SC326143**

ZNF780A (NM_001142577) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF780A (NM_001142577) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF780A
Synonyms:	ZNF780
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326143 representing NM_001142577.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGTCCATGGATCAGTGACATTCAGGGATGTGGCCATTGACTTCTCTCAGGAGGAGTGGGAGTGCCTG
CAGCCTGATCAGAGGACCTTGACAGGGATGTGATGTTGGAGAACTACAGCCACCTGATCTCACTGGCA
GGAAAGTTCCATTTCTAAACCAGATGTAATTACGTTACTAGAGCAAGAGAAAGAGCCCTGGATGGTTGTA
AGGAAAGAAACAAGCAGACGGTATCCAGATTTGGAGTTAAATATGGACCTGAGAAAGTATCTCCAGAA
AATGATACCTCTGAAGTAAATTTACCCAAACAGGTTATAAAGCAAATAAGTACAACCTCTTGGCATTGAG
GCCTTTTATTTAGAAATGACTCAGAATATAGACAATTTAGGGACTACAGGGATATCAAGAAGGAAAT
ATCAATCAAAGATGATCAGCTATGAAAACTGCCTACTCATACTCCTCATGCTTCTTATTTGCAAT
ACACATAAACCGTATGAATGTAAGGAATGTGGGAAATACTTTAGTCGTAGTGCAAATCTTATTCAGCAT
CAGAGTATTCATACTGGAGAGAAACCTTTGAATGTAAGGAGTGTGGGAAAGCCTTCGACTTCACATA
CAATTTACTCGACATCAGAAATTTCACTGCTGAGAAACCTTTGAATGTAACGAATGTGGAAAGGCC
TTTAGTCTTCTTACCCTGCTTAATCGCCATAAGAACATTCACACAGGTGAGAAATGTTTGAATGTAAG
GAATGTGGGAAGTCCTTTAATCGTAGCTCAAACCTTGTTCAACATCAGAGTATTCATTCTGGTGTAAAA
CCATATGAATGTAAGGAGTGTGGGAAAGGCTTTAATCGTGGTGCACACCTTATTCAGCATCAGAAATTT
CATTCCAATGAGAAACCTTTGTATGTAAGGAATGTGGGATGGCCTTTCGATATCATTACCAACTTATT
GAACATTGCCAAATTCATACTGGTGAGAAACCTTTGAATGTAAGAATGTGGAAAGGCGTTTACTCTT
CTGACAAAGCTTGTTCGACATCAGAAGATTCATACTGGTGAGAAACCTTTGAATGCAGGGAATGTGGG
AAGGCCTTTAGTCTTCTCAACCAGCTTAATCGCCATAAGAACATTCACACAGGTGAAAAACGCTTTGAA
TGTAAGGAATGTGGGAAGTCCTTTAATCGTAGCTCAAACCTTGTTCAACATCAGAGTATTCATGCTGGT
ATAAAACCATATGAATGTAAGGAGTGTGGGAAAGGCTTTAATCGTGGTGCACACCTTATTCAGCATCAG
AAAATTCATTCCAATGAGAAACCTTTTGTATGTAGGGAATGTGAGATGGCCTTTAGATATCATTGCCAA
CTTATTGAACATTCTCGAATTCATACTGGTGACAAGCCATTTGAATGTCAAGACTGTGGGAAAGCCTTC
AATCGTGGCTCAAGCCTTGTTCAACATCAGAGTATTCACACTGGTGAGAAAGCCCTATGAATGTAAGGAG
TGTGGGAAGGCTTTTAGACTTTACCTACAACCTTTCCCAACATCAGAAACTCACACAGGTGAAAAACCA
TTTGAATGTAAGGAATGTGGGAAATTTTCGTCGTGGTTCAAATCTTAATCAACATCGAAGTATTCAT
ACTGGAAGAAACCTTTGAATGTAAGGAATGTGGGAAAGCCTTTGACTTCATATGCACCTTATTCGA
CATCAGAAATTCGATACTGGTGAGAAACCTTTGAATGTAAGGAGTGTGGGAAAGCCTTTGACTTCAT
ATGCAACTTATTCGACATCAGAAATTCGATACTGGTGAGAAACCTTTGAATGTAAGGAATGTGGAAAG
GTTTTTAGTCTTCCACCCAGCTTAATCGCCATAAGAACATTCACACAGGTGAGAAGGCATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001142577

Insert Size: 1929 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001142577.1](#)

RefSeq Size: 3615 bp

RefSeq ORF: 1929 bp

Locus ID: 284323

Cytogenetics: 19q13.2

MW: 74.6 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
 Transcript Variant: This variant (1) encodes the longest isoform (a).