

Product datasheet for **SC326154**

ZNF534 (NM_001143938) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF534 (NM_001143938) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF534
Synonyms:	KRBO3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326154 representing NM_001143938.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCTTACTCAGGGCAATTGTCATTACGCGATGTGGCCATAGAATTCTCTCAGGAGGAGTGGAAA
TGCCCTGGACCCTGGGCAGAAAGCTTTATACAGGGACGTGATGTTAGAGAACTACAGGAACCTGGTCTCC
CTAGGAATCTGTCTTCTGACCTGAGTGTACCTCCATGTTGGAGCAAAAGAGAGATCCCTGGACTCTG
CAGAGTGAAGTAAAAATAATAACAATCCAGATGGCAGGGAGTGCATCAAAGGTGTGAACACAGAGAAG
AGTTCTAAATTGGGAAGCAGTGCAGAAACAAGTCACTTAAAAATCAACATGGATTAACCTTTCAGTTA
CATCTGACTGAATGGCAGCCATTTCAAGCTGTAAGGAATATTTATGGATGAAGCATGTTGAGAAATCT
ATCAGTGACAATCTTCAGTTTCACCAGTTCAAAATAAGTTTTTTCAGTGTCAAACCCATATTTTTAAT
AACTACAGAAATGATTTTCTTTTTCTACACTCCACAAGAACAGAAAGTACACATTAGGGAAAAG
CCTTATGGATGTAATGAGCATGGGAAAGTCTTCAGAGTGTCTCAAGCCTTACTAACCGTCAAGTAATC
CACATTGCAGATAAACTTACAAATGTAGTGACTGTGGCGAGATCTTTAGTAGCAATTCAACTTTGCA
CAACATCAACGAATCCATACTGGAGAGAAGCCTTACAAATATAATGAATGTGGCAAAGTCTTCAATCAG
AATTCACACCTTGCAACAATCAGAAAATTCATACTGGACAGAAAACCTTACAATAACAAAGAATGTGGG
AAAGTCTTTAGTACCATGCCTACCTTGACAGCATAGGAAAATTCATACTGGAGAGAAGCCTTACAAA
TGTAGTGAATGTGGCAAAGCATTAGTGTGTGTTCCAGTCTTACTGCTCATCTTGTAAATCCATACTGGA
GAGAAACCTTATGATTGTAAGGAATGTGGCAAGTCTTCAGGCATAAGTCTTCCCTAACCACTCATCAG
ACAGTTCATACTGGAGAGAGACCATACAAATGTAATGAATGTGGCAAGGGTTTAGTCGAATTCATTTC
CTTGAAGGCATCGGAAAGTTCATACTGGAGAGAAAACCTTACAAATGTAATGAGTGTGGCAAGTCTTT
ATTGGCAATTCAGCCTTGACGACATAGGAAAATTCATACTGGGGGGAGGCGTTACAAATGTAATGAA
TGTGGCAAAGCATTTAGAACGTGTTAGATCTCACTGCCCATCTTCTAATCCATACTGGAGAGAAAACCT
TACGAATGTATAGACTGTGGCAAGTCTTCAGGCACAAGTCTTCCCTAACCTATCACTGTAGAATTCAT
ACTGGAGAGAAGCCTTACAAATGTAACGAATGTGGCAAGTCTTCAGTCAGAATTCAAACCTTCAACGA
CATAGGAAAATTCATACTGGAGAGAAGCTTTACAAATGTAATGAATGTGGCAAGTCTTCCGTCAGAAT
TCACACCTTGCAACAATAGGGATTCATACTGGAGAGAAGCCTTACAGTTGTAATGAATGTGGCAAG
GTCTTCCGTCGGAATTCACACCTTGTGCGACATAGGAATGTTCACTGGAGAAAAGCCTTACAGTTGT
AATGAATGTGGCAAGTCTTCAGTCGGAATTCACACCTTGCAGCATAGGAATTCATACTGGAGAG
AAGCCTCACAGTTGTAATGAATGTGGCAAGTCTTCAGTCGGAATTCACACCTTGCAGCATAGGAAA
ATTCATACTGGAGAGAAGCTTTACAAATGTAATGAATGTAGCAAGTCTTCAGTCGGAATTCACGCCTT
GCACAACATAGGAATTCATACTGGAGTGAAGCCTTACAGTTGTAATGAATGTGGCAAGTCTTTAGT
AAAAATTCATCCTAGTACAACATTGCAGTATTCATACCAGAGAGAAGCCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001143938

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001143938.3</u>
RefSeq Size:	5203 bp
RefSeq ORF:	1986 bp
Locus ID:	147658
UniProt ID:	<u>Q76KX8</u>
Cytogenetics:	19q13.41
MW:	75.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.