

Product datasheet for **SC326155**

Epsin 1 (EPN1) (NM_001130071) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Epsin 1 (EPN1) (NM_001130071) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPN1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC326155 representing NM_001130071.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGTGATCAGAGCTGGCTCTGGAACCAAGGCTGCTCCAGGAGTAAGGAGCCAGTCTTTGCATGCAGT
GTTGAAAAGGTAAATGTCCCTCTTGCTGAGCGAGCACCTGGCACACAGCAGGGACCCAGGCAGTGGG
GCTGTTAGGTTCCCTATCTCTCCTGAGCCTTGGGCTTCCGCTATCCTTGGGACGTCTGGTCTTGGCT
TCCCCGGTTCTTCTGCCGCCCTGGATGCGGTGACCTGCCAGCACCTGCCGCAGCCTTCGTCCGGGAGT
CGCCCCATCTCTCCACGCATCGGGGCCCTGTGCCCTTGTGCTGCAGCCGGGCACCATGTCGACCTCG
TCCTTGAGGCGCCAGATGAAGAACATCGTCCAACTACTCAGAGGCGGAGATCAAGGTTTCGAGAGGCC
ACGAGCAATGACCCCTGGGGCCCATCCAGCTCCCTCATGTCAGAGATTGCCGACCTCACCTACAACGTT
GTCGCCTTCTCGGAGATCATGAGCATGATCTGGAAGCGGCTCAATGACCATGGCAAGAACTGGCGTCAC
GTTTACAAGGCCATGACGCTGATGGAGTACCTCATCAAGACCGGCTCGGAGCGGTGTGCGAGCAGTGC
AAGGAGAACATGTACGCCGTGCAGACGCTGAAGGACTTCCAGTACGTGGACCGCGACGGCAAGGACCAG
GGCGTGAACGTGCGTGAGAAAGTAAGCAGCTGGTGGCCCTGCTGCGCGACGAGGACCGGCTGCGGGAA
GAGCGGGCGCACGCGCTCAAGACCAAGGAAAAGCTGGCACAGACCGCCACGGCCTCATCAGCAGCTGTG
GGCTCAGGCCCCCTCCCGAGGCGGAGCAGGCGTGGCCGAGAGCAGCGGGGAGGAGGAGCTGCAGCTC
CAGCTGGCCCTGGCCATGAGCAAGGAGGAGGCCGACCAGGAGGAGCGGATCCGTCGCGGGGATGACCTG
CGGCTGCAGATGGCAATCGAGGAGAGCAAGAGGGAGACTGGGGGCAAGGAGGAGTCTGCCCTCATGGAC
CTTGCTGACGCTTTCACGGCCCCAGCTCCTGCCCGACCACAGACCCCTGGGGGGGCCAGCACCCATG
GCTGCTGCCGTCCCACGGCTGCCCCACCTCGGACCCCTGGGGGCCCCCTGTCCCTCCAGCTGCT
GATCCCTGGGGAGGTCCAGCCCCACGCCGGCTCTGGGGACCCCTGGAGGCTGCTGCCCTGCAGGA
CCCTCAGTTGACCCTTGGGGTGGGACCCAGCCCTGCAGCTGGGGAGGGGCCACGCGCTGATCCATGG
GGAAGTCCGATGGTGGGTCCCGTCACTGGGCCCTCAGCCTCCGATCCCTGGACACCGGCCCCGGCC
TTCTCAGATCCCTGGGAGGGTCACTGCCAAGCCAGCACCAATGGCACAAACAGCAGCCGGGGATTTC
GACACGGAGCCCGACGAGTTCTCTGACTTTGACCGACTCCGCACGGCACTGCCGACCTCCGGGAGCAGC
GCAGGAGAGCTGGAGCTGCTGGCAGGAGAGGTGCCGGCCGAAGCCCTGGGGCGTTTGACATGAGTGGG
GTCAGGGGATCTCTGGCTGAGGCTGTGGGCAGCCCCCACCTGCAGCCACACCAACTCCCACGCCCCC
ACCCGGAAGACGCGGAGTCAATCCTGGGGCCCAATGCAGCCCTCGTCGACCTGGACTCGCTGGTGGC
CGCCCGGGCCCCACGCCCTGGAGCAAGGCCTCCAACCCCTTCTGCCAGGCGGAGGCCAGCCACT
GGCCCTTCCGTCACCAACCCCTTCCAGCCCGCCCTCCCGGACGCTCACCTGAACCAGCTCCGCTC
AGTCCGTGTCCTCCCGTCCCTGGAGCGCCACCCACGTACATCTCTCCCTTGGCGGGGGCCCTGGCCTG
CCCCCATGATGCCCCCGGGCCCCCGGCCCAACTAATCCCTTCTCCTATATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001130071

Insert Size: 1989 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130071.1](#)

RefSeq Size: 2964 bp

RefSeq ORF: 1989 bp

Locus ID: 29924

UniProt ID: [Q9Y6I3](#)

Cytogenetics: 19q13.42

Protein Pathways: Endocytosis

MW: 69 kDa

Gene Summary: This gene encodes a member of the epsin protein family. The encoded protein binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Loss of function of this gene is associated with reduced tumor growth and progression in certain cancer types. [provided by RefSeq, Mar 2016]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).