

Product datasheet for **SC326159**

ZNF237 (ZMYM5) (NM_001142684) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF237 (ZMYM5) (NM_001142684) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF237
Synonyms:	HSPC050; MYM; ZNF198L1; ZNF237
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326159 representing NM_001142684.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAAAAATGTTTCAGTGGGAGGATTAGAGTTGACTGAACAGACTCCTGCTTTATTAGGGAATATGGCC
ATGGCAACTAGTCTCATGGACATAGGGGATTCATTTGGTCATCCAGCTTGTCCTTTAGTCAGTAGATCT
AGGAACCTCACCAGTGGAAGATGATGATGATGATGATGATGTTGTGTTTATTGAATCTATACAACCTCCT
TCAATTTCTGCTCCAGCAATAGCTGATCAAGAACTTCATATTTGCATCATCAAAAAATGAAAAGCCT
CAAGGAAATTATTCTGTAATTCCTCTTCAAGAGATTTGGCATCTCAGAAAGGAAATATAAGTGAG
ACAATTGTTATTGATGATGAAGAGGACATAGAAACAAATGGAGGAGCAGAGAAAAAGTCTTCTGTTTT
ATCGAATGGGGACTTCTGGAATAAAACAAACCAACGATTTGGATTTCTCCACTCCAGTCTTTCA
AGAAGTAAGACCAAGACTGGAGTAAGACCTTTTAAACCTGGTAGAATGAATGTGGCAGGAGACTTATTT
CAGAATGGAGAATTTGCAACTCATCATAGTCTGATTCTTGGATCTCCAGTCAGCTTCATTCCCAGT
AATCAGAAACAGCCAGGGGTGGATTCTTATCACCAGTGGCCTTACTTCGTAAGCAGAATTTCCAGCCT
ACAGCCCAACAACAACCTACTAAACAGCTAAATCACTTGTGCAAAATGCAAAAAGCCTTTACAGAAG
GGACAGACAGCTTATCAACGAAAAGGATCAGCTCACCTCTTTGTCTACCACTGCCTTTCTCTCTTC
TCTCATAAACGTAATAAAACACACGAAGCATAATATGTAAGAAAGATGCATCTACAAAGAGGCTAAT
GTCATTCTTCCAGTAGAATCAAGCAATCCTTCCAAGATTTTATAGTACATCTTGTGTTGTCTCCCTGT
GAAAACAACCTGGAATCTTAAAAAAGGAGTTTTTAATAAGTCAAGATGTACAATTTGTAGTAAATTAGCA
GAGATTCGCCATGAAGTCAGCGTAAATAATGTAACACATAAACTGTGCAGTAACCATTTGCTTTAATAAG
TACAGATTGGCCAATGGTCTAATAATGAAGTCTGTGAACACTGTGGAGAGTACATGCCTAGTAAGAGT
ACTGAAAACAACATCCTGGTGTGAGGTGAGGAGAGAGATTTGCTGCCAAAGTTGTATTAACGAA
TATAAACAGATGATGGAACAAATCAAAAAATTAACAGCATCAGAAATAGAAAAAGGAATGCTTTT
AGAGAAGAAAATGAGAAACAATTATATGGATCGTCAAATACACTTTTAAAAAATAGAGGGAATCCCA
GAAAAAAGGAAAAGACTTCACAGCTACAGCTTTCAGTAGAATGTGGCACGGATACATTACTGATCCAA
GAGAATGTGAATTTACCTCCTTCTTCCAGTCAACCATAGCTGATACATTTCAAGAGCAACTGGAAGAG
AAAAATTTTGAAGACTCCATTGTACCAGTTGTGCTTTCCGCAGATCCAGGTACGTGGCCCGAATTTTG
AATATTAACAACGGGACACTCTTGTGAAAATGTTCCGCCTCAAGTAAGAAATTTTAACTTTCAAAAA
GATAATACAGGGAGGAAGTTTTAGAACTTATTATACAGGATCTTCCAAATGGTAAAAAACCACT
AGATCCTGGTTACTTTTCAACCTCAAAAGATTCTGTGTTTGTCTATATTGCAAACTCTTTGGGGAA
GGAAAAATCAACTGAAAAATGAAAATGGGTGCAAGATTGGCAACATTTATCACATATTCTTAGTAAA
CATGAAGAAAGTGAATGCACGTCAACAATAGTGTAAAGTACTCAAAATTAATCTGATCTGAAAAAA
AATAAAGCTATTGATGCTGCAGAACACAGATTATATGAAAATGAGAAAAATGATGGTGTGCTGTTGTTG
TACACGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001142684

Insert Size: 2010 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142684.1</u>
RefSeq Size:	3268 bp
RefSeq ORF:	2010 bp
Locus ID:	9205
UniProt ID:	<u>Q9UJ78</u>
Cytogenetics:	13q12.11
MW:	74.8 kDa
Gene Summary:	Functions as a transcriptional regulator.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (3). There are no full-length transcripts representing the exon combination of this variant in human; it is predicted based on partial transcript alignments, on alignment of the macaque mRNA AB179330.1, and on homology support for the encoded protein. Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The extent of this RefSeq transcript is supported by transcript alignments.