

## Product datasheet for **SC326162**

### EPB4IL2 (EPB41L2) (NM\_001135554) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	EPB4IL2 (EPB41L2) (NM_001135554) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPB41L2
Synonyms:	4.1-G; 4.1G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC326162 representing NM\_001135554.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACTACTGAAGTAGGCTCTGTGTCTGAAGTGAAGAAGGACTCTAGCCAGTTAGGAACAGATGCAACC
AAGGAAAAACCTAAAGAAGTAGCAGAAAAATCAGCAGAATCAGTCTCCGATCCAGAGGAGGAAAAAGGT
TCCCAGCCACCTCTGCAGCTGAAAGCCAAAGTAGTCTACGCCGCCAGAAGAGAGAGAAGGAAACATCG
GAGAGCAGGGTATTTCTCGTTTCATACCGCCATGGCTTAAGAAGCAAAAAGTCATATACCTTAGTAGTG
GCCAAAGATGGAGGAGATAAAAAAGAGCCTACCCAAGCTGTTGTTGAAGAACAGGTCTTAGATAAAGAG
GAACCCCTTCCAGAAGAACAGAGACAGGCTAAGGGTGATGCTGAAGAAATGGCTCAGAAGAAACAAGAG
ATTAAGTTGAAGTCAAGGAAGAAAAACCTCAGTGAGCAAGGAAGAAAAACCTCAGTGAGCAAAGTG
GAGATGCAGCCTACTGAATTAGTAAGTAAGGAGAGAGAAGAGAAGGTAAAAGAAACACAGGAAGACAAA
TTAGAAGGAGGAGCAGCAAAAAGGGAGACCAAGGAAGTGCAGACCAATGAGCTGAAAGCAGAGAAGGCA
TCTCAAAAAGTCACCAAGAAGACCAAACTGTCCAGTGTAAAGTGACCCTTTAGATGGCACCGAATAC
AGCTGTGACCTGGAGAAACATGCCAAGGACAAGTGTTATTTGACAAAGTGTGTGAACACCTCAATCTC
TTGGAGAAAGACTACTTTGGACTTTTGTTCAGAAAAGCCCTGAGCAGAAAAACTGGTTAGATCCTGCT
AAAGAAATAAAGAGACAACCTGAGAAACCTTCCATGGCTATTCACCTTTTAAATGTGAAGTTTATCCTCCT
GATCCTTCTCAATTGACTGAAGATATCACCAGATACTTCTGTGCCTTCAGCTCCGGCAGGACATTGCC
TCTGGCCGCTGCCCTGCTCTTTTGTGACTCATGCTCTCCTGGGATCCTACACCTGCAGGCTGAACTT
GGTGACTATGACCAGAAGAACATGGCAGCATCGACCTCAGTGAATTCAGTTTGGCCCTACTCAGACT
AAGGAGCTGGAAGAGAAGGTGGCAGAGCTGCACAAAACCCACAGGGGCTTATCGCCAGCACAAAGCTGAT
TCCAGTTCTTAGAAAAATGCAAAGAGGCTTCCATGTATGGTGTGACCTACATCATGCCAAGGACTCA
GAAGGTGTGGACATCAAGCTGGGCGTGTGTGCTAATGGACTTCTCATTTACAAAGACAGACTGCGAATC
AATCGTTTTGCTTGGCCGAAAATCTTAAAAATTCCTATAAACCGCAGTAACCTTACATTAAGTACAGA
CCGGCAGAGCTGGAACAGTTTGAGAGTACCATTGGATTCAAACCTGCCAAACCACCGGGCAGCGAAAAAGA
CTATGGAAGTGTGCGTGGAGCATCATACTTTCTACAGGCTTGTTCCTCAGAGCAGCCACCAAAAGCC
AAGTTCTGACCTTGGGGTCCAATTTGCTATAGTGGCCGCACCAAGCACAGACCCGCCAGGCCAGC
ACCCTCATAGATAGGCCAGCACCACTTTGAGCGCACTTCTAGTAAACGGGTCTCCAGGAGTCTAGAT
GGAGCTCCGATTGGTGCATGGACCAAAGTCTTATGAAGGATTTTCTGGCGCTGCTGGGAGATTTCA
GCCTATGGACCTGGACTTGTGAGCATTGCCGTGGTACAAGATGGGACCGCAGGAGGGAAGTGAAGAGC
CCAACTAAAGCCACATTTGCAGCTCATTGAAGGAAAGACTGTAAAAGGTGGAATTTCTGAAACAAGA
ATTGAGAAACGCATTGTGATCACAGGAGATGGAGATTTGATCATGACCAGGCACTGGCTCAGGCGATC
AGGGAAGCCAGAGAGCAGCACCCCTGACATGTCCGTACAAGAGTGGTGGTACACAAGAAACAGAGTTG
GCTGAGGAAGGGGAAGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001135554

**Insert Size:** 2022 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001135554.1](#)

**RefSeq Size:** 3473 bp

**RefSeq ORF:** 2022 bp

**Locus ID:** 2037

**UniProt ID:** [O43491](#)

**Cytogenetics:** 6q23.1-q23.2

<b>Protein Pathways:</b>	Tight junction
<b>MW:</b>	76.1 kDa
<b>Gene Summary:</b>	Required for dynein-dynactin complex and NUMA1 recruitment at the mitotic cell cortex during anaphase (PubMed:23870127).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) lacks multiple consecutive in-frame exons in the CDS, compared to variant 1. The resulting isoform (b) lacks an internal segment, compared to isoform a. Variants 2 and 3 encode the same isoform (b).