

Product datasheet for **SC326173**

AGAP5 (NM_001144000) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGAP5 (NM_001144000) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGAP5
Synonyms:	CTGLF2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC326173 representing NM_001144000.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGGGAACATACTGACCTGTTGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAACAGGGGTCG
GTGTGTCCCTCTGAATCTGAGATCTATGAGGCAGGAGCTGGGACAGGATGGCAGGAGCGCCCATGGCT
GCTGCTGTGCAGCCTGCTGAGGTGACCGTTGAAGTTGGTGAGGACCTCCACATGCACCACATTTCGTGAC
CAGGAGATGCCTGAAGCTTTGGAGTTTAGCCTTTCTGCCAATCCAGAGGCAAGCACAATATTCCAGAGG
AACTCTCAAACAGATGCTTTGGAGTTTAAACCTTCTGCCAATCCAGAGGCAAGCACAATATTCCAGAGG
AACTCTCAAACAGATGTTGTAGAAATAAGAAGAAGCAACTGTACAAACCATGTATCTACTGAGCGTTTC
AGTCAACAATACAGCTCGTGTTCGACAATATTCCTTGATGACAGCACAGCCAGCCAGCATTATCTTACA
ATGACAATAATATCTGTGACCTTGAGATACCTCATCATATCACACAAAGAGATGCAGATAGATCTTTG
AGCATACCTGATGAACAGTTACACTCATTTCGGGTTTCCACCGTGCACATTACGAAGAACAGAAATGGA
GGTGGGAGTTTAAATAACTATTCTCTCCATTCCATCGACTCCAGCACCAGCCAGGAGGACCCTCAG
TTCAGTGTTCCTCCACTGCCAACACACCCACCCCGTTTGAAGCGGTCCATGCCGTGGTCCAACCTG
TTTACATCTGAGAAAGGGAGTCACCCAGACAAAGAGAGGAAAGCCCCGGAGAATCATGCTGACACCATC
GGGAGCGGCAGAGCCATCCCCATTAAACAGGGCATGCTCTTAAAGCGAAGTGGGAAATGGCTGAAGACA
TGGAAAAAGAAATACGTACCCCTGTGTTCCAATGGCGTGCTCACCTATTATTCAAGCTTAGGTGATTAT
ATGAAGAATATTCATAAAAAAGAGATTGACCTTCGGACATCTACCATCAAAGTCCAGGAAAGTGGCCA
TCCCTAGCCACATCGGCTGTGCACCCATCTCCAGCTCTAAAAGCAATGGCCTATCCAAGGACATGGAC
ACCGGGCTGGGTGACTCCATATGCTTCAGCCCCAGTATCTCCAGCACCACCAGCCCCAAGCTCAACCCG
CCCCCTCTCCTCATGCCAATAAAAAGAAACACCTAAAGAAGAAAAGACCAACAACCTTTATGATTGTG
TCTGCCACTGGCCAAACATGGCACTTTGAAGCCACGACGTATGAGGAGCGGGATGCCTGGGTCCAAGCC
ATCCAGAGCCAGATCCTGGCCAGCTGCAGTCATGCGAGAGCAGTAAAAGCAAGTCCAGCTGACCAGC
CAGAGCGAGGCCATGGCCCTGCAGTCGATCCAAAACATGCGTGGGAACGCCCACTGTGTGGACTATGAG
ACCCAGAATCCTAAGTGGGCCAGTTTGAAGTGGGAGTCTCATGTATTGAATGCTCAGGAATCCAC
CGCAGTCTTGGCACCCGCTTTCCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGG
AAGTTATGTATCTATTGGCAATGACCTAGCCAACAGCATCTGGGAAGGGAGCAGCCAGGGGCAGACA
AAACCCTCAGTAAAGTCCAGAGGGAAGAGAAGGAACGGTGGATCCGTTCCAAATATGAGGAGAAGCTC
TTTCTGCCCCACTACCCTGCACTGAGCTGTCCCTGGGCCAGCACCTGCTGCGGGCCACCGCTGATGAG
GACCTGCAGACAGCCATCCTGCTGCTGGCACATGGCTCCCGTGAGGAGGTGAACGAGACCTGTGGGAG
GGAGACGGCTGCACGGCGCTCCATCTGGCCTGCCGAAGGGGAATGTGGTCTGGCACAGCTCCTGATC
TGGTACGGGGTGGACGTATGGCCCGAGATGCCACGGGAACACAGCGCTGACCTACGCCCGGCAGGCC
TCCAGCCAGGAGTGCATCAACGTGCTTCTGCAGTACGGCTGCCCGCAGGAGTGCGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001144000

Insert Size: 2061 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001144000.1</u>
RefSeq Size:	2426 bp
RefSeq ORF:	2061 bp
Locus ID:	729092
UniProt ID:	<u>A6NIR3</u>
Cytogenetics:	10q22.2
MW:	75.5 kDa
Gene Summary:	Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]