

## Product datasheet for **SC326178**

### Invadolysin (LMLN) (NM\_001136049) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Invadolysin (LMLN) (NM_001136049) Human Untagged Clone
Tag:	Tag Free
Symbol:	LMLN
Synonyms:	GP63; INV; IX14; LMNL1; MSP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC326178 representing NM\_001136049.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGTAACGACGCTCGGCCGAAGATGGCGCCGAATGGGGCGGAGGAGTGGTTACTCGGGCTCAGGC
CCGGCCCGGAGCCGGTGGCGCTGGAGCGGTCTGTGGGTCCGAAGCGTTTACTCCTGTTGGCGGG
CTCCGGCCAGCGCCACATCTACTCCCGTCTCCTTGGGAGTCCCTCCCTGCCGGCACCACGTCCCC
TCTGACACTGAGGTCATAAATAAAGTTCATCTTAAGGCAATCATGTGGTCAAGAGAGATGTTGATGAG
CATTAAAGAATCAAGACTGTCTATGATAAAAGTGTGAAGAGTTGCTCCCTGAGAAAAAGAATCTTGTA
AAGAACAAGCTTTTCCACAAGCGATTTCTTATTTAGAGAAGACTTTTCAGGTCCGTCGACCTGCCGGC
ACTATCTTACTTAGCAGACAATGTCAACAAACCAATACCTCCGGAAGGAAAACGATCCTCACAGGTAC
TGCACCGGGGAGTGTCCGCACACACAAAGTCCGCCCCGTTATTGTTCTGAGGAACATCTCCAGCAA
TGCCGGGTCTACCGTGGGGTAAGTGGCCTCATGGAGCAGTGGGTGTCCAGACCAAGAAGGCATCTCA
GATGCAGACTTTGTTCTTACGTTGGTGTCTGGCCACCGAGAGATGCAGCCATGAAAACATCATCTCT
TATGCAGCCTATTGTGACAGGAAGCAAACATGGACAGGCCAATAGCAGGATATGCTAACCTGTGTCCA
AATATGATCTCTACCCAGCCTCAGGAGTTTGTGGGATGCTGTCCACAGTGAAAACATGAGGTTATTCAT
GCCCTGGGTTTCTGCTGGGCTGTTTGCATTCTACCATGATAAAGATGGAAATCCTCTCACTTCAAGA
TTTGCAGATGGCCTCCCACCTTTAATTATAGTCTGGGATTATATCAATGGAGTGATAAAGTAGTTTCA
AAAGTGGAGAGATTATGGGATGTTTCAGATAAATAAGATAGTTCGTACACTGTGTATCTCCTGGTAAAG
CCTCGTGTGTTGAGGAAGCAGAAAACATTTGATTGTCAGTTCAGAGGGAATGGAAGTGAAGTGAAGT
CAAGTGGTGTGGGCACTGAGCTCAACCATTTGGGAAAAAGGTTATTAGAGAATGAAGCGATGACTGGT
TCTCACACTCAGAATCGAGTACTCTCGAATCACTCTGGCATTAAATGGAGGACACTGGCTGGTATAAA
GAAAATTACAGCATGGCTGAGAAGTTAGACTGGGGCCGAGGAATGGGCTGTGACTTTGTGAGGAAGAGC
TGTAATTTCTGGATTGATCAGCAGAGACAAAAGAGACAGATGCTGAGCCCTTACTGTGACACGCTCAGA
AGTAACCCACTGCAGCTAACTTGCAGACAGGACCAGAGAGCAGTTGCCGTGTGTAATTTGCAGAAGTTC
CCTAAGCCTTTACCACAGGAATACCAGTACTTTGATGAACTCAGTGAATACCTGCAGAAGATTTGCCT
TATTATGGTGGCTCCGTGGAATTTGCTGACTACTGCCCTTTCAGTCAGGAATTCAGTTGGCATTAAAGT
GGTGAATATCAGCGCAGCTCAGATTGTAGAATATTGGAAAATCAACCAGAAAATTTTAAAGAACTATGGC
GCTGAAAAGTATGGACCTCATTCCGTTTGTCTAATTCAGAAAATCAGCATTGTTATGGAGAAGTGTGAG
AGGAAGCTGAGTTACCCAGACTGGGAAGCGGATGCTATCAGGTTTCTTGTCTCCTCAAGTCTGAAA
GTTTGGTCCAAGATACTTCATATTTGTGTAGTCGGGCTGGGCAGGTCCCTCCCTGTGAGTATCCAGATG
AATGGCTGGATTACGATGGAACCTGCTCTGCCATCATGTTGGGACTTCTGTGAGCTCTGTCCCTCCA
GAAACAGATCCTCCAGCCACTAACCTGACCCGAGCTCTGCCACTTGATCTTTGTTCTGTTCTCGAGC
CTGGTGGTACCCTCTGGCTTCTGCTAGGCAATCTGTTTCTCTGCTGGCTGGATTTCTTCTGTGTATA
TGGCACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001136049

**Insert Size:** 2079 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001136049.2](#)

**RefSeq Size:** 7154 bp

**RefSeq ORF:** 2079 bp

**Locus ID:** 89782

**UniProt ID:** [Q96KR4](#)

**Cytogenetics:** 3q29

**Protein Families:** Druggable Genome, Protease, Transmembrane

**MW:** 78.1 kDa

**Gene Summary:** This gene encodes a zinc-metallopeptidase. The encoded protein may play a role in cell migration and invasion. Studies of a similar protein in *Drosophila* indicate a potential role in mitotic progression. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.